

US Patent & Trademark Office

Patent Public Search | Text View

United States Patent Application Publication

20250255828

Kind Code

A1

Publication Date

August 14, 2025

Inventor(s)

PRADO; ISAIAS et al.

LYOPHILIZED LIPID NANOPARTICLES AND METHODS OF THEIR USE

Abstract

The invention is directed to the field of therapeutic formulations, in particular to lyophilization of a therapeutic cargo molecule, such as RNA. The invention provides a method for lyophilization of a molecule. The present disclosure further describes a lyophilized composition obtainable by the inventive method, a pharmaceutical composition, a vaccine, a therapeutic and a kit or kit of parts. Moreover, the disclosure herein provides a novel lyophilization excipient that protects the composition from degrading when, for example, lyophilizing RNA. The use of the inventive method further includes the manufacture of a composition that can be used after lyophilization with equivalent therapeutic effect and composition integrity.

Inventors: PRADO; ISAIAS (Novato, CA), BHAMBHANI; AKHILESH (Novato, CA), CHENG; BENJAMIN (Novato, CA), ZHANG; HAIRUI (Novato, CA), WANG; XIAOWEI (Hillsborough, CA), GARRIPELLI; VIVEK KUMAR (Novato, CA)

Applicant: Ultragenyx Pharmaceutical Inc. (Novato, CA)

Family ID: 96661324

Assignee: Ultragenyx Pharmaceutical Inc. (Novato, CA)

Appl. No.: 19/047781

Filed: February 07, 2025

Related U.S. Application Data

us-provisional-application US 63551205 20240208

Publication Classification

Int. Cl.: A61K9/51 (20060101); A61K9/19 (20060101); A61K31/7105 (20060101); A61K38/17 (20060101); A61K38/45 (20060101); A61K38/47 (20060101)

U.S. Cl.:

CPC A61K9/5123 (20130101); A61K9/19 (20130101); A61K31/7105 (20130101); A61K38/1709 (20130101); A61K38/45 (20130101); A61K38/47 (20130101);

Background/Summary

CROSS-REFERENCE TO RELATED APPLICATIONS [0001] This application claims the benefit of and priority to U.S. Provisional Application No. 63/551,205, filed 8 Feb. 2024, the entire disclosure of which is hereby incorporated herein by reference for all purposes.

BACKGROUND OF THE INVENTION

[0002] Transfection of nucleic acid therapeutics and other agents has been accomplished by encapsulating the active molecules in lipid nanoparticles. Drawbacks of this methodology include the inability to store compositions for later use because of degradation of the nanoparticles or the encapsulated contents therein. For example, compositions of lipid nanoparticles that encapsulate nucleic acid molecules may be stable for only a few minutes or hours at 25° C., and only a few days or weeks at 4° C. Further drawbacks include the need for very low temperature storage of the lipid nanoparticle compositions.

[0003] One way to provide for long-term storage of a therapeutic composition is to prepare a lyophilized form, which can be stored and reconstituted to provide a formulation for administration of the therapeutics.

[0004] However, it has not been possible in general to generate lyophilized forms of lipid nanoparticles containing nucleic acid agents, so that the lipid nanoparticle can be regenerated with the nucleic acid agent encapsulated to form a stable formulation. The lyophilization process can destroy the nanoparticles and/or the nucleic acid agents. Some methods have involved chemically attaching protective groups or components to the lipid nanoparticles, or to the nucleic acid agent, which is disadvantageous. Other methods may use liposomes as an alternative to LNP, without providing for encapsulation and/or stability of the nucleic acid agents.

[0005] There is a continuing need for compositions and methods to provide lyophilized forms of nanoparticles that can be reconstituted with favorable properties, including transgene activity or potency, particle size, storage stability, and shelf-life to deliver various nucleic acid agents.

[0006] What is needed are methods, compositions and compounds for forming stable solutions or suspensions of lipid nanoparticles that can be stored in solid lyophilized forms, where the nanoparticles encapsulate a cargo molecule or composition.

[0007] In gene therapy as well as in many other therapeutically relevant biochemical and biotechnological applications such as enzyme replacement therapies, nucleic acids are used for therapeutic and diagnostic purposes. As an example, rapid progress has occurred in recent years in the field of gene therapy and promising results have been achieved. Nucleic acids are therefore regarded as important tools for gene therapy and prophylactic and therapeutic vaccination against, for instance, genetic disorders as well as infectious and malignant diseases. Gene therapy has also shown to be an effective therapeutic against uncommon diseases that affect a small portion of the population (rare and ultra-rare diseases).

[0008] Nucleic acids, both DNA and RNA as well as small antisense oligonucleotides (ASOs), have been used widely in gene therapy, either in naked or in complexed form. In this context, the application of nucleic acids and particularly of RNA for therapeutic treatment of genetic disorders

is becoming ever-increasingly important.

[0009] The application of RNA represents a favored tool in modern molecular medicine, which also exhibits some superior properties over the application of DNA. As generally known, transfection of DNA molecules may lead to serious complications. For example, application of DNA molecules bears the risk that the DNA integrates into the host genome. Integration of foreign DNA into the host genome can have an influence on the expression of host genes and can trigger the expression of an oncogene or the inactivation of a tumor suppressor gene. Furthermore, an essential gene—and, as a consequence, the product of such an essential gene—may also be inactivated by the integration of the foreign DNA into the coding region of the gene.

[0010] The result of such an event may be particularly dangerous if the DNA is integrated into a gene, which is involved in regulation of cell growth. Notwithstanding the risks associated with its application, DNA still represents an important tool. However, these risks do not occur if RNA, particularly mRNA, is used instead of DNA. An advantage of using RNA rather than DNA is that no virus-derived promoter element has to be administered in vivo and no integration into the genome may occur. Furthermore, the RNA, in order to exert its function, does not need to overcome the barrier to the nucleus. However, a main disadvantage of the use of RNA is its instability. For instance, RNA is susceptible to hydrolysis by ubiquitous ribonucleases or by divalent cations and is typically rapidly degraded, e.g., after a few hours or days in solution. Rapid degradation occurs even in the absence of RNases, e.g., when RNA is stored in solution at room temperature for a few hours or days. To avoid such rapid degradation, RNA (in solution) is typically stored at -20°C . or even -80°C . and under RNase free conditions to prevent degradation of the RNA. Such storage conditions, however, do not sufficiently prevent a loss of function over time. Additionally, applying such conditions is very cost-intensive, especially for shipping and storage, e.g., whenever such low temperatures have to be guaranteed. Thus, there is a need for stable biological therapeutics using, for example, siRNA, mRNA, or other RNA-based modalities.

[0011] One further method for stabilization of RNA comprises lyophilization or freeze-drying of the RNA. Lyophilization is a method known and recognized in the art, which is used to enhance storage stability of temperature sensitive biomolecules. During lyophilization, a solvent, such as water, is typically removed from a frozen sample via sublimation. However, a nominal amount of water may remain.

[0012] The process of lyophilization is usually characterized by a primary and a secondary drying step. During the primary drying step, free, i.e., unbound, water surrounding the biomolecule and optionally further components, sublimates from the frozen solution. Subsequently, water, which is bound by the biomolecule on a molecular basis, may be removed in a secondary drying step by adding thermal energy. In both cases, the hydration sphere around the biomolecule is lost.

[0013] During lyophilization, a sample containing a biomolecule is initially cooled below the freezing point of the solution and accordingly of the water contained therein. As a result, the water freezes. Depending, amongst other parameters, on temperature, cooling rate (freezing rate), and the time for freezing, crystals may be formed. This exerts physical stress on the biomolecule and other components of the solution, which may lead to damage of the biomolecule such as—in the case of a nucleic acid—breakage of strands, structural alteration, etc. Furthermore, due to the decrease of volume and loss of the hydration sphere, autocatalytic degradation processes are favored, e.g., by traces of transition metals. In addition, the concentration of traces of acids and bases can result in significant changes of the pH value. Lyophilization involves two types of stress, namely freezing and drying. Both types of stress are known to damage nucleic acid molecules.

[0014] In order to reduce or minimize damage to the molecule a number of cryoprotectants are discussed for lyophilization purposes to prevent these damages. In this context, cryoprotectants are understood as excipients, which allow influencing the temperature characteristics of the mixture. Cryoprotectants can also partially or totally replace the hydration sphere around a molecule and

may thus at least partially prevent catalytic and hydrolytic processes. Lyophilization may increase the stability of DNA under long-term storage, but may also cause some damage due to the initial lyophilization process, potentially through changes in the DNA secondary structure, breaks of the nucleic acid chain(s) or the concentration of reactive elements such as contaminating metals.

[0015] One solution for improving the storage capabilities of lipid nanoparticle formulations is to manufacture the lipid nanoparticle formulations as a lyophilized product that can be subsequently reconstituted prior to administration. Lyophilized lipid nanoparticle compositions can be stored at more practical temperatures, allowing more convenient modes of distribution and administration.

[0016] Although lyophilization technologies have existed for several decades, the application of these technologies does not translate well to lipid nanoparticle formulations, which lose several of their desired properties including low polydispersity, small particle size, high percentage of encapsulation, and in vivo efficacy upon reconstitution after conventional lyophilization techniques. Therefore, improved formulations and methods are needed to allow for initial lyophilization processes to maintain structural integrity of the encapsulated cargo molecule and also maintain therapeutic effectiveness upon reconstitution.

BRIEF SUMMARY OF THE INVENTION

[0017] The present disclosure provides lyophilized compositions and methods of making lyophilized compositions.

[0018] The disclosure is directed to the field of therapeutic formulations, in particular to lyophilization of a therapeutic cargo molecule, such as nucleic acid agents (e.g., mRNA). The disclosure provides methods for lyophilization of a cargo molecule encapsulated inside a lipid nanoparticle (LNP). The present disclosure further describes lyophilized compositions obtainable by the disclosed methods, such as pharmaceutical compositions, vaccines, and other therapeutics. Moreover, the disclosure herein provides a novel approach to produce mRNA/LNP solution pre-lyophilization that protects the composition from degrading when, for example, lyophilizing RNA. The use of the disclosed methods further includes the manufacture of a composition that can be used after lyophilization with equivalent therapeutic effect and composition integrity.

Description

BRIEF DESCRIPTION OF DRAWINGS

[0019] FIG. 1 is a graphical illustration of intensity weighted particle distribution.

[0020] FIG. 2 is a graph showing Small-angle X-ray Scattering (SAXS) patterns of the measurement samples.

[0021] FIG. 3 is a cartoon depiction of proposed sample of a lipid nanoparticle containing mRNA.

[0022] FIG. 4A is a process flow chart for preparing a sample for lyophilization.

[0023] FIG. 4B is a graphical illustration of size and PDI under various buffer formulations.

[0024] FIG. 5 is a flow chart showing an mRNA lyophilization production protocol.

[0025] FIG. 6 is a graphical illustration of lyophilization cycle profile for Round 1 screening completion of primary drying process.

[0026] FIG. 7 is a graphical illustration of lyophilization cycle profile for Round 2 screening completion of primary drying process.

[0027] FIG. 8 is a graphical illustration of lyophilization cycle profile for Round 3 screening completion of primary drying process.

[0028] FIG. 9 shows supplemental data of Round 1 samples pre-lyophilization.

[0029] FIG. 10 shows supplemental data of Round 1 samples post-lyophilization.

[0030] FIG. 11 shows supplemental data of Round 1 samples reconstituted with Milli-Q water.

[0031] FIG. 12 shows supplemental data of Round 2 samples pre-lyophilization.

[0032] FIG. 13 shows supplemental data of Round 2 samples post-lyophilization.

[0033] FIG. **14** shows supplemental data of Round 2 samples reconstituted with Milli-Q water.
[0034] FIG. **15** shows supplemental data of Round 2 samples reconstituted with TNS Buffer.
[0035] FIG. **16** shows supplemental data of Round 3 samples pre-lyophilization.
[0036] FIG. **17** shows supplemental data of Round 3 samples post-lyophilization.
[0037] FIG. **18** shows supplemental data of Round 3 samples reconstituted with Milli-Q water.
[0038] FIG. **19** shows size stability results observed for lots 230523A and B stored 1 week at 40° C., and 6 months at 25° C., 5° C. and -20° C.
[0039] FIG. **20** shows stability study results for lots 230523A and B by polydispersity index [PDI].
[0040] FIG. **21** shows stability results by encapsulation efficiency observed for lots 230607A and B after 1 week at 40° C., 1 weeks at 25° C., and 6 months at 5° C. and -20° C.
[0041] FIG. **22** shows stability results by encapsulation efficiency for lots 230523A and B.
[0042] FIG. **23** shows purity study results for lots 230607A and B after 2 weeks at 40° C. and 25° C., and 6 months at 5° C. and -20° C.
[0043] FIG. **24** shows purity study results for lots 230523A and B after 7 days at 40° C., and 6 months at 25° C., 5° C., and -20° C.

DETAILED DESCRIPTION OF THE INVENTION

[0044] The present disclosure provides lyophilization methods that result in a preservation of lipid nanoparticle integrity, integrity of the encapsulated molecule, the particle size of the lipid nanoparticles within an acceptable degree of pre-lyophilized particle size, with an improved excipient formulation. The methods stem from the discovery that specific combinations of excipients can be added to a pretreated suspension of the mRNA and eventually to nanoparticles prior to subjecting the suspension to a lyophilization process. In addition, lyophilization parameters are employed in combination with these excipients to achieve high quality lyophilized lipid nanoparticle product. The lyophilized product is easily reconstituted and readily administered as a pharmaceutical.

Definitions

[0045] “Cargo molecule” means peptides or proteins, DNA, RNA, carbohydrates, lipids or chemically devised molecules of natural or non-natural origin, such as metals, dyes, and nanoparticles that can be encapsulated and/or delivered as a carrier for a therapeutic agent and/or a therapeutic agent. Non-limiting examples of cargo molecules include DNA, RNA, miRNA, mRNA, siRNA, synergistic activation mediator (SAM), circular RNA, ASO, or a combination thereof.

[0046] The term DLS as used herein means Dynamic Light Scattering.

[0047] “Effective” refers to capability of a therapeutic or diagnostic agent, or an amount thereof (e.g., mass, volume, dosage, concentration, and/or time period) to achieve one or more desired result(s).

[0048] “Encapsulation efficiency” (EE) as used herein is the fraction of initial drug that is encapsulated by the nanoparticles (NPs).

[0049] “Excipient” as used herein refers to an organic or inorganic ingredient, natural or synthetic inactive ingredient in a formulation, with which one or more active ingredients are combined.

[0050] “Nanoparticle” as used herein refers to a lipid containing particle that can be formed using a T-mixing approach, microfluidic approach, jet-mixing, solvent emulsion, spray drying, or precipitation in bulk or microfluids, wherein the solvent is removed to no more than an insignificant residue, leaving an electron dense particle of size ≤ 200 nm, preferably ≤ 150 nm, ≤ 100 nm (which may, or may not, be hollow or have a liquid filled interior such as a liposome) polymeric particle, unlike a micelle whose form is dependent upon being present in an aqueous solution.

[0051] “N/P: Molar ratio” is defined as the ratio of nitrogen on ionizable lipid to Phosphate on the nucleic acid backbone.

[0052] “Prevention” or “preventing” means to administer a composition to a subject or a system at

risk for or having a predisposition for one or more symptom caused by a disease or disorder to cause cessation of a particular symptom of the disease or disorder, a reduction or prevention of one or more symptoms of the disease or disorder, a reduction in the severity of the disease or disorder, the complete ablation of the disease or disorder, stabilization or delay of the development or progression of the disease or disorder.

[0053] The term “protein” “polypeptide” or “peptide” refers to a natural or synthetic molecule comprising two or more amino acids linked by the carboxyl group of one amino acid to the alpha amino group of another.

[0054] The term R.sub.h as used herein means Hydrodynamic Radius.

[0055] The term R.sub.g as used herein means Radius of Gyration.

[0056] The term SAM refers to synergistic activation mediator that is an engineered protein complex for the transcriptional activation of endogenous genes.

[0057] The term SAXS as used herein means Small-angle X-ray Scattering.

[0058] The present disclosure provides improved lyophilized compositions and methods of making lyophilized compositions. Specifically, the present disclosure is directed to the field of molecule formulations, in particular lyophilization of a molecule, such as a nucleic acid agent. The disclosure provides methods for lyophilization of a molecule. The present disclosure further describes lyophilized compositions obtainable by the disclosed methods, pharmaceutical compositions, vaccines, and other therapeutics. Moreover, the disclosure herein provides a novel lyophilization excipient that protects the composition from degrading when, for example, undergoing the process of lyophilization. The use of the disclosed methods further includes the manufacture of a composition that can be used after lyophilization with equivalent therapeutic effect and composition integrity.

[0059] As described in further detail herein, the method of making a composition comprises a LNP and a cargo molecule encapsulated within the LNP. In exemplary embodiments, the composition is capable of being lyophilized. The method includes combining an aqueous solution of the cargo molecule and an ionizable lipid composition contained in an organic solvent solution, thereby forming a lipid/cargo molecule suspension.

[0060] In example embodiments, the cargo molecule that is combined with the lipid composition comprises at least one nucleic acid, such as DNA or RNA. In related embodiments, the cargo molecule comprises micro RNA (miRNA), messenger RNA (mRNA), short interfering RNA (siRNA), short hairpin RNA (shRNA), antisense oligonucleotide (ASO), synergistic activation mediator (SAM), circular RNA, or a combination thereof.

[0061] The lipid composition can include an ionizable lipid composition used to form an ionizable organic solution. In these embodiments the ionizable lipid composition comprises an ionizable cationic lipid, such as a tertiary amine that is capable of deprotonation under neutral conditions and is positively charged in pH conditions below the acid-dissociation constant (pK.sub.a) of the lipid. In other embodiments, the ionizable organic solvent solution comprises at least one of an ionizable lipid, DSPC, cholesterol, PEG-DMG and combinations thereof. In some related embodiments, the ionizable organic solvent solution comprises a N/P molar ratio of ≥ 2.5 , ≥ 6.0 , or ≥ 9.0 . It should be appreciated that the combining the aqueous solution of the cargo molecule and a lipid composition with an ionizable organic solvent solution forms a lipid/cargo molecule suspension.

[0062] In example embodiments, excipients are added to the lipid/cargo molecule suspension. In other related embodiments, the method of making the composition comprises excipients where the excipients are added to the aqueous solution of the cargo molecule prior to being combined with the lipid composition. In at least these example embodiments, the excipients are selected from a group comprising citrate, tris, arginine, NaCl, sucrose, PVA, phosphate, HEPES, trehalose, KCl, acetate, bis-tris, histidine, glucose, lactose, raffinose, alanine, asparagine, proline, glutamic acid, methionine, threonine, and combinations thereof. In other example embodiments, excipient(s) comprise at least one positively charged and at least one negatively charged compound or an amino

acid combination. The excipients that can be used within the scope of the invention, include but are not limited to at least two amino acids selected from Ala, Arg, His, Asp, Pro, Glu, Met, Thr, and combinations thereof. In other embodiments, the excipient comprises a hydrophilic polymeric compound. In other related embodiments, the excipient comprises a plasticizer.

[0063] In example embodiments, the aqueous solution and/or lipid composition comprises an organic solvent. Examples of organic solvents within the scope of the present disclosure include at least one of ethanol, DMSO, acetone or other alcohol containing DSPC, cholesterol, PEG-DMG, an ionizable lipid, and combinations thereof. In the methods of making the lipid/cargo molecule suspension, the process comprises diluting the lipid/cargo molecule suspension that is combined with the organic solvent to lower the concentration of organic solvent to $\leq 10\%$ of the total composition. In alternative embodiments involving the methods of making the lipid/cargo molecule suspension, the process comprises diluting the lipid/cargo molecule suspension that is combined with the organic solvent to lower the concentration of organic solvent to $\leq 3\%$ of the total composition.

[0064] The dilution step comprises a first dilution step to increase the pH to about 6 and a second dilution step to increase the pH to between 6.8 and 7.8. The first dilution step comprises at least one of phosphate, HEPES, tris, bis-tris, acetate, citrate, NaCl, KCl, glutamic acid, arginine, histidine, methionine, glucose, sucrose, lactose, trehalose, raffinose, and combinations thereof. The second dilution buffer comprises at least one buffer compound selected from amino acids, sugars, salts and combinations thereof. These dilutions steps may further comprise the step of adding a stabilizer to the aqueous solution of the cargo molecule.

[0065] The stabilizer is added prior to mixing the cargo molecule with the organic solvent. The stabilizer is optionally present in any one of the at least one dilution steps.

[0066] When forming the lipid/cargo molecule suspension, the pH of the composition is increased to at or above 5.0. When modifying the pH, the adjustment to at or above a pH of 5.0 can be performed before the dilution step (b). In other alternative embodiments, when modifying the pH, the adjustment to at or above a pH of 5.0 can be performed after the dilution step (b). In either pH modification step being performed prior to, or after the dilution step, the pH can be adjusted to at least 6. In other embodiments, the pH is adjusted to between 6.9 and 7.9.

[0067] As described above, once the lipid/cargo molecule suspension is formed, the lipid/cargo molecule suspension can undergo the process of lyophilization. In these embodiments lyophilizing a composition comprises the lipid nanoparticle having the cargo molecule encapsulated therein, thereby forming a lyophilized cargo molecule/LNP, the method comprising freezing the lipid/cargo molecule suspension. Once the freezing step is completed, a further drying step can be performed on the frozen lipid/cargo molecule suspension. During the drying process, the drying step can be performed by a process selected from a group comprising spin-freeze drying, continuous freeze-drying, sublimation, desorption, vacuum foam drying, microwave vacuum drying, spray drying, evaporation and combinations thereof. During the drying step, the target sublimation is achieved when the moisture is present at an amount of $\leq 5\%$.

[0068] Once lyophilization is completed, the cargo molecule/LNP is tested for potency and/or purity yield. Potency and/or purity yield of the lyophilized composition is $\geq 70\%$, $\geq 80\%$, $\geq 90\%$, or $\geq 95\%$ when compared to the cargo molecule/LNP prior to lyophilization. Other characteristics of the lyophilized cargo molecule/LNP include the determination of particle size and/or encapsulation efficiency. Particle size undergoing the process of lyophilization described herein, includes; < 100 nm, < 90 nm, < 80 nm, or < 75 nm, while encapsulation efficiency of the lyophilized cargo molecule/LNP is about $\geq 80\%$, or about $\geq 90\%$.

[0069] In some exemplary embodiments, the cargo molecule within the cargo molecule/LNP comprises a therapeutic agent at a starting concentration of ≤ 2 mg/ml under acidic conditions with a pH under 5, preferably between pH 3-4.5. In some exemplary embodiments, the cargo molecule within the cargo molecule/LNP comprises mRNA as the therapeutic agent at a starting

concentration of ≤ 1 mg/ml under acidic conditions with a pH under 5, preferably between pH 3-4.5. In some exemplary embodiments, the cargo molecule within the cargo molecule/LNP comprises mRNA as the therapeutic agent at a starting concentration of < 0.55 mg/ml under acidic conditions with a pH under 5, preferably between pH 3 and pH 4.5. In some exemplary embodiments, the cargo molecule within the cargo molecule/LNP comprises mRNA as the therapeutic agent at a starting concentration of ≤ 0.22 mg/ml under acidic conditions with a pH under 5, preferably between pH 3 to pH 4.5. In an example embodiment, a cargo molecule/LNP comprises mRNA encoding amylo-alpha-1, 6-glucosidase, 4-alpha-glucanotransferase (AGL), wherein the mRNA encoding AGL is codon-optimized. In some embodiments, the AGL is human AGL (i.e. hAGL). Example cargo molecules that can be utilized within the scope of the present disclosure can be found in U.S. Pat. No. 11,377,643, the contents of which are incorporated herein in its entirety. These cargo molecules can be used for ameliorating, preventing, delaying onset, or treating a disease or condition associated with reduced activity of amylo-alpha-1, 6-glucosidase, 4-alpha-glucanotransferase (AGL) in a subject. Other cargo molecules can be used to treat diseases, such as progressive familial intrahepatic cholestasis (PFIC), including PFIC-1, PFIC-2 and PFIC-3. For example, a cargo molecule/LNP may comprise a mRNA encoding human bile salt export pump (human BSEP). The mRNA may be an ABCB11 mRNA encoding human BSEP. The mRNA may be codon-optimized.

[0070] In the described process of forming the cargo molecule/LNP compositions, the cargo molecule, in a non-lyophilized state, prior to encapsulation, is dissolved in a buffer under acidic pH conditions below the pKa of the lipid composition. The buffer may optionally comprise a salt between 25 mM and 100 mM. In some embodiments, the salt is NaCl. The buffer may optionally comprise 50 mM NaCl.

[0071] The process of lyophilization of the present disclosure provides an improvement in cargo molecule/LNP integrity, such that the step of reconstitution of the lyophilized cargo molecule/LNP composition can be performed with a diluent.

[0072] The lipid nanoparticle in the composition after being lyophilized and then reconstituted is increased in size between 3-18% or to a size $\leq 18\%$ relative to the size of the lipid nanoparticle prior to lyophilization. In other embodiments, the lipid nanoparticle in the composition after being lyophilized and then reconstituted is increased in size no more than 18% relative to the size of the lipid nanoparticle prior to lyophilization.

[0073] Prior to lyophilization, the defined cargo molecule, where the cargo molecule is RNA, has an encapsulation efficiency (EE) of the RNA of at least 70%, 80%, 85%, or 90%. Importantly, the encapsulation efficiency (EE) of the composition after being reconstituted has an average variance of no more than 20% of the lyophilized or aqueous form of the composition. In other embodiments, the encapsulation efficiency (EE) of the composition after being reconstituted has an average variance of between 10-20% of the lyophilized or aqueous form of the composition.

[0074] The methods of forming a lyophilized composition can also include the process of preparing a lyophilized composition comprising lipid nanoparticles with reduced percentage of water in a lipid nanoparticle (LNP) core, the method comprising: [0075] (a) combining an aqueous solution of a cargo molecule with a lipid composition, thereby forming a lipid/cargo molecule suspension; [0076] (b) mixing the lipid/cargo molecule suspension with an organic solvent; [0077] (c) adding an excipient, wherein the excipient is added before or after the aqueous solution of the cargo molecule is combined with the lipid composition; [0078] (d) diluting the lipid/cargo molecule suspension having the organic solvent to lower the concentration of organic solvent to $\leq 10\%$ of the total composition by volume; thereby forming a lipid nanoparticle composition comprising lipid nanoparticles comprising a shell of the lipid composition and a core comprising the cargo molecule in aqueous solution with the organic solvent; [0079] (e) lyophilizing the lipid nanoparticle composition; and [0080] (f) drying the lipid nanoparticle composition, thereby forming the lyophilized composition with lipid nanoparticles having reduced percentage of water in the lipid

nanoparticle core.

[0081] In some embodiments, the process of reducing water in the LNP core as described herein can reduce the moisture content of the lyophilized composition to $\leq 5\%$ by weight. Moreover, when performing step (f), the drying step can be performed by a process selected from spin-freeze drying, continuous freeze-drying, sublimation, desorption, vacuum foam drying, microwave vacuum drying, spray drying, evaporation, and combinations thereof. The drying step of at least this embodiment occurs at a temperature between 30°C . to -50°C . The drying step is generally performed at a time range of between 10-3000 minutes.

[0082] The process may further comprise reconstituting the lyophilized composition in a diluent. The diluent, in some embodiments, includes saline, water, 50-500 mM NaCl, arginine, D5W, sucrose, polyvinyl alcohol (PVA), tris, or any combinations thereof.

[0083] Samples (see Table 1) were prepared using the prescribed methods and were analyzed using small-angle X-ray Scattering (SAXS) measurement on three different samples including 1 control (220316) and 2 samples (210825) & (1119A). Buffer was also measured for background subtraction. SAXS results suggests that 220316 forms Lipid Nanoparticles (LNPs) with random lipid distribution in the core and with PEG lipids on the surface of LNP. Following the introducing of mRNA, 1119A and 210825 show bilayer structure, which represents the interaction between the lipids and mRNA. Bilayer stacks were also observed in SAXS measurements which is consistent with TEM images (not shown).

TABLE-US-00001 TABLE 1 Samples received

Sample Number	Sample ID #	Sample Description
1	220316	Control
2	1119A	Sample
3	210825	Sample
4	220808	Buffer

[0084] Sample preparation is performed by DLS measurement using a Zetasizer Ultra Red system with a side detector mode (90°). SAXS measurements were performed using an Empyrean Nano system with a ScatterX78 SAXS accessory with a PIXCel 3D detector. Approximately 40 μL of sample were loaded in a 1.0 mm (O.D.) quartz capillary tube. The covered scattering vector, ($q=4\pi/\lambda \sin \theta/2$ where θ is the scattering angle and λ is the incident X-ray wavelength= 1.5418 \AA) is from $\sim 0.006\text{ \AA}^{-1}$ to 0.5 \AA^{-1} . Sample scattering and buffer scattering were measured in the same capillary tube for proper background subtraction. EasySAXS software was used for background subtraction which utilized the attenuated incident beam intensity as the normalization factor. SAXS data analysis is performed on SASView 5.0.5.

[0085] For sample 220316, control, a combined Core-shell model with a polymer chain model was used to fit the SAXS pattern. For samples 210825 and 1119A, a combined core 2 shell model with a Lorentzian peak. A combined 2 core shell model, as used herein is described in *Small* 18, no. 9 (2022) by Cui et al, entitled “Mechanistic studies of an automated lipid nanoparticle reveal critical pharmaceutical properties associated with enhanced mRNA functional delivery in vitro and in vivo”.

DLS Measurement

[0086] Referring now to FIG. 1, DLS was performed to verify the sample freshness after shipping. The intensity weighted particle size distribution results are shown. The Z-averaged hydrodynamic radius ($R_{\text{sub.h}}$) as well as polydispersity index (PDI) are listed in Table 2. The $R_{\text{sub.h}}$ measured by Malvern Panalytical (MP) was within the required specification. Therefore, the SAXS measurements were performed.

TABLE-US-00002 TABLE 2 Hydrodynamic radius and PDI for each sample

Sample	$R_{\text{sub.H}}$ (NM)	$R_{\text{sub.H}}$ UX (NM)	PDI MP	PDI UX
220316	52.5 ± 1.2	38	0.19 ± 0.01	0.139
210825	48.7 ± 1.4	40.0	0.15 ± 0.02	0.125
1119A	40.5 ± 0.1	36.5	0.08 ± 0.03	0.08

SAXS Measurement

[0087] SAXS was performed and the results are shown in FIG. 2. The fitting parameters are listed in Table 3. For all the samples measured, the size is comparable and sample 1119A and 210825 has the largest and smallest PDI, respectively. For sample 220316, since no mRNA is added, no lipid bilayer (phosphate-alkyl chain-phosphate) is formed. The core of the 220316 is a lipid mixture

while the PEG-lipid was primarily located at the surface with PEG faced to the solvent/buffer forming a PEG shell. The proposed LNP structure of sample 220316 is shown in FIG. 3(a). The structural parameters were also listed where radius of the core is ~ 19 nm and the radius of gyration of PEG is ~ 1 nm. Once lipids are mixed with mRNAs, i.e. sample 210825 and 1119A, the bilayer structure is formed.

[0088] In addition, these bilayers further form stacks which include water layers as shown in FIG. 3(b) where the small circles represent phosphate group with alkyl chain and water layer as yellow and blue area. Here, one stack is defined as phosphate-alkyl chain-phosphate-water. Therefore, in FIG. 3(b) the proposed structure has 1 stack plus one phosphate. The ratio between diffraction peak and the core-shell structure reveals the number of bilayer stacks. The electron density in the bilayer stack is higher in 210825. Moreover, the electron density is larger in the bilayer stack than 1119A. Based on the RiboGreen assay, results were indicative that the electron density in the bilayer stack is higher in 210825.

[0089] SAXS results confirmed the bilayer structure formed with presence of mRNA. In addition, the structural parameters can be extracted from the SAXS pattern by fitting the data. Moreover, SAXS results are consistent with the TEM images and the RiboGreen assay. Referring still to FIG. 3 (Proposed structure from Small Angle X-ray Scattering (SAXS)), analysis suggests the presence of bilayer stacks with water in between them. The small circles represent phosphate group with alkyl chain and water layer illustrated as yellow and blue area. While not wanting to be bound by any particular theory, a skilled artisan, in view of the detailed description herein, could determine that replacing water with stabilizers during the mRNA/LNP process early on would be essential to enable successful lyophilization, e.g., excipients that serve as a water replacer (e.g. sucrose and trehalose) and/or known for charge shielding (e.g. NaCl) and Arg/Glu salt pairs. The desired ranges were then determined empirically, and it was unexpectedly found that the presence of Glu/Arg, in presence or absence of His, performed better than Arginine alone or in absence of any amino acid. Moreover, the packing density defect during the dehydration process manifested as an increase in particle size upon reconstitution. Thus, minimizing shock (due to charge and/or osmotic imbalances) during reconstitution may enable retention of particle size post-lyophilization. In prior lyophilization protocols, in-use stability for non-lyophilized formulations (i.e. frozen liquid) revealed that diluting out stabilizer resulted in dramatic change in size and loss in recovery post reconstitution.

TABLE-US-00003 TABLE 3 Fitting parameters used for SAXS analysis

Core	Shell 1	Shell 2	Peak	Scale	Sample	(NM)	(NM)	(Å)	sup.-1	(NM)	Peak	220316	19.3	1.0	1.1	N/A	N/A	N/A	N/A	(Fixed)	210825	22.0	2.4	N/A	1.0	0.115	5.5	0.56	(Fixed)	(Fixed)	1119A	22.6	2.4	N/A	1.0	0.115	5.5	0.42	(Fixed)	(Fixed)
Radius	Thickness	Polymer Thickness	Position	Spacing	CS/Scale																																			

[0090] The methods described herein provide an improved method of successfully lyophilizing mRNA/LNP compositions with retention of mRNA purity, relative potency, % EE and particle size with low polydispersity index (PDI) as measured using Capillary Electrophoresis (CE), In-cell Western (ICW), RiboGreen Assay and Dynamic Light Scattering (DLS).

[0091] It should be appreciated that limited shelf-life of mRNA/LNP is driven in part by presence of water in the LNP core that results in, for example, hydrolysis of phosphodiester bond in the mRNA backbone and lyophilization serves as an appealing method to improve the shelf-life of mRNA/LNP. Additionally, LNPs are self-assembled multi-complex systems and stresses induced during the freeze-drying (such as freezing, vacuum, dehydration, etc.) processes are known to impact the physicochemical properties (i.e. particle size, PDI, EE % etc.) and ultimately mRNA-LNP functionality (i.e. potency or efficacy). Moreover, the addition of stabilizers early during the process, starting with mRNA pre-mix, and optionally optimizing the composition may serve as a means to overcome the challenges associated with dehydration of mRNA/LNP using freeze-drying. Stabilizers such as sucrose and trehalose serve as water replacers; however, their use in mRNA pre-

mix has not been observed before during mRNA/LNP process. Similarly, the presence of additional stabilizer that can form salt bridges such as Glu+Arg as part of mRNA dilution buffer have not been observed before. Surprisingly, the presence of stabilizer such as sucrose, glutamate and arginine in mRNA dilution buffer favorably reduced the particle size to <70 nm (versus the original process with a particle size of ~80 nm).

[0092] As shown in FIGS. 4A and 4B (see also Table 4), two different approaches to LNP production prior to lyophilization are shown. In the traditional approach (referred to as pC-1 in Table 4), mRNA/LNP are fabricated and lyoprotectants are added prior to lyophilization. In the alternate approach (referred to as pC-2 in Table 4), selected cryoprotectants and lyoprotectants are introduced in the mRNA dilution buffer itself. The bar graph (FIG. 4B) demonstrates that presence of Glu+Arg+Sucrose provides lowest particle size prior to lyophilization while keeping PDI <0.1. TABLE-US-00004

TABLE 4 Approaches to LNP production	Step	pC-1	pC-2
mRNA Dilution Buffer	6.5 mM Citrate, 50 mM 5 mM Glutamic Acid, 5 mM NaCl, pH 4.0	Histidine, 25 mM Arginine, 15% Sucrose, pH 4.0	Dilutions Buffer 1 45 mM Phosphate, pH 6.0 50 mM Histidine, 9% Sucrose, pH 6.0 Dilution Buffer 2 100 mM Tris, 15% Sucrose, pH 7.6
Diafiltration Buffer	50 mM Tris, 15% Sucrose, pH 7.4	Concentration Excipient and 50 mM Tris, 500 mM Arginine, 250 mM NaCl, 15% Adjustment Buffer 1 Concentration Sucrose, 0.5% PVA, pH 7.4	Adjustment Concentration 50 mM Tris, 100 mM Arginine, 50 mM NaCl, 15% Adjustment Sucrose, 0.1% PVA, pH 7.4

[0093] While not reciting any single factor as being critical, several factors are required to minimize the impact of freeze-drying on mRNA/LNP physicochemical (such as size, PDI, EE % etc.) and functional properties (based on % purity, potency, etc.) such as: (a) Process of making mRNA/LNP by ensuring presence of stabilizers during (mixing, dilution step 1, dilution step 2, lyophilization once mRNA/LNPs are formed, and reconstitution); (b) Approach to mixing with microfluidics yielding more consistent particles than T-mixer; (c) Starting mRNA concentration wherein lower starting concentration of mRNA results in lower particle size pre-lyophilization; (d) Presence of salt wherein [NaCl]≥25 mM is needed maintain particle <80 nm; (e) Presence of Glu+Arg+sucrose yields particle size <80 nm; (f) Lyophilization process parameters also impacted particle size and corresponding potency and purity results; (g) Presence of Tris and Sucrose by itself is not sufficient to keep particle size <100 nm and presence of Arg+PVA10+NaCl is needed with earlier addition of excipients (pre-mixing and/or dilution 1: PC2), preferred over late addition (after dilution step 2; PC1); and (h) Custom reconstitution media such as formulation buffer, saline, half saline performed better than water alone in retaining particle size.

[0094] Stabilization of mRNA/LNP remains a challenge, as exemplified by cold-chain storage requirement of COVID vaccines. Liquid formulations are preferred from a customer convenience perspective; however, they pose sterility assurance and stability challenges often requiring cold-chain storage. Lyophilization serves to overcome the stability challenges and in early-stage clinical program provides an approach for early entry into clinic by attaining stable product with flexibility on dosing and drug concentration by using custom reconstitution media and volume. Thus, there remains a consensus in the field that freezing and lyophilization are the most common approach to retain long-term storage and shelf life of lipid nanoparticles. As discussed above, the stresses induced by lyophilization often results in a dried mRNA/LNP of suboptimal characteristics. Under certain conditions lyophilization of mRNA/LNP results in a suboptimal particle size (often ≥100 nm) with PDI, often ≥0.1 or a loss of EE (in some cases ≤90%) or a compromise on potency and/or purity. Comparative lyophilization processes, such as those found in US 2022/0047519 documents the use of specialized excipient to retain potency. Similarly, the addition of 30% ethanol during the reconstitution step results in a size comparable to fresh LNPs, however, that process is neither convenient nor practical for use in the clinic (see Ball et al, Int J Nanomedicine 2016 December 30:12:305-315.doi: 10.2147/IJN.S123062. Use of trehalose and sucrose seems to alleviate the problem; however, the particle size was high (~200 nm) with low EE (60-70%). Similarly, previous

studies (see Muramatsu, Mol Ther. 2022 May 4; 30 (5): 1941-1951)) demonstrated a particle size of ~96-99 nm for lyophilized and reconstituted Luc mRNA/LNPs, while the latter reference demonstrates particle size of lyophilized LNP to be greater than 150 nm.

[0095] The present disclosure exemplifies the process of making stable lyophilized mRNA/LNP formulations and their resulting composition of matter using Test samples as cargo molecule/LNP composition. Test sample is a lipid encapsulated therapeutic mRNA expressing a gene of interest (GOI) formulated at 1 mg/mL in 50 mM HEPES, 50 mM NaCl, 9% Sucrose, 5% Glycerol pH 7.4 and stored as frozen liquid stable formulation. Two different approaches are shown using multiple compositions. In the first approach (referred to as pC-1), mRNA/LNP are fabricated and lyoprotectants are added prior to lyophilization. In the alternate approach (referred to as pC-2), selected cryo- and lyoprotectants are introduced in the mRNA dilution buffer itself.

TABLE-US-00005 TABLE 5 pC-1 and pC-2 mRNA-LNP lyophilization Size PDI Encapsulation DP Pre Post Pre Post Pre Post Formulation Lyo Lyo Δ Lyo Lyo Δ Lyo Lyo Δ Potency pC-1 50 mM Tris, 83 95 11 0.135 0.111 -0.024 98% 95% -3% NA 100 mM Arg, 0.1% PVA10, 50 mM NaCl, 9% Sucrose, pH 7.4 20 mM Tris, 83 97 14 0.110 0.130 0.020 98% 96% -2% NA 100 mM Arg, 0.1% PCA10, 50 mM NaCl, 9% Sucrose, pH 7.4 50 mM Tris, 83 92 9 0.116 0.131 0.015 98% 94% -4% NA 100 mM Arg, 0.1% PVA10, 25 mM NaCl, 9% Sucrose, pH 7.4 pC-2 50 mM Tris, 71 101 30 0.037 0.101 0.064 99% 86% -13% 222% 9% Sucrose, pH 7.4 50 mM Tris, 71 95 24 0.049 0.117 0.068 99% 87% -12% 191% 15% Sucrose, pH 7.4 50 mM Tris, 71 83 12 0.034 0.068 0.034 99% 95% -4% 164% 100 mM Arg, 0.1% PVA10, 50 mM NaCl, 9% Sucrose, pH 7.4 50 mM Tris, 71 82 11 0.031 0.065 0.034 99% 94% -5% 144% 100 mM Arg, 0.1% PVA10, 50 mM NaCl, 15% Sucrose, pH 7.4

[0096] mRNA/LNP size pre-lyophilization is lower when excipients are introduced early in mRNA dilution buffer (pC-2) versus (pC-1). For the similar process pC-2, the presence of Arginine, PVA10 and NaCl enables greater retention of particle size upon lyophilization. When final formulations were optimized, the formulation included 50 mM Tris, 100 mM Arginine, 0.1% PVA10, 50 mM NaCl, 9% Sucrose, and a pH of 7.4 for both the process condition (i.e. pC-1 and pC-2, respectively).

[0097] In alternative embodiments, mRNA dilution buffer components, such as salt, can be altered to further improve the invention described herein. For example, effects of [NaCl] on Particle Characteristics: mRNA diluted buffer ranged from 5-75 mM NaCl. After 30 min. incubation period between Dilution 1 Pool and Dilution 2 Pool, the LNPs were evaluated for size, PDI, and encapsulation. PDI is lowest when [NaCl] is between 35-50 mM NaCl (See Table 6 below). As shown, the effect of salt concentration of mRNA/LNP particle size, PDI and encapsulation efficiency demonstrates that [NaCl]≥25 mM results in acceptable mRNA/LNP particle size.

TABLE-US-00006 TABLE 6 mRNA dilution buffer components [NaCl] in mRNA Encapsulation Dilution Buffer (mM) Size (nm) PDI Efficiency (%) 75 81 0.088 97 62.5 78 0.121 98 50 77 0.069 97 35 73 0.038 96 25 79 0.156 99 15 96 0.190 98 5 140 0.184 95

[0098] In other related embodiments, modifications of Amino acids: mRNA/LNP size follows the order Glu/Arg (+His)<Arg alone<no amino acid can also be used to provide improvements in LNP/cargo molecule compositions. In other related embodiments, sugars, such as sucrose, can be used to show improvement in size and PDI.

[0099] In related embodiments, the modification of an alcohol can also be used to modify the EE. As shown in Table 7, the percentage of ethanol during mixing is shown. Particle characteristics when % EtOH ranges from 10-30% are stable (5% Ethanol condition has slightly lower encapsulation). As % EtOH decreases, mRNA must be decreased to maintain the N: P ratio. Thus, mRNA/LNP are stable when % EtOH ranges from 10-30% while 5% EtOH condition results in a slightly lower encapsulation efficiency.

TABLE-US-00007 TABLE 7 Modification of alcohol content Encapsulation Hydration Time % EtOH Size (nm) PDI Efficiency (%) (1 g Batch) 30 81 0.108 99 6.7 min 25 77 0.069 97 8.1 min 20

75 0.062 97 10.1 min 15 72 0.049 95 13.5 min 10 79 0.035 92 20.2 min 5 79 0.053 89 40.4 min [0100] As shown, lower starting concentration results in lower particle size pre-lyophilization (e.g. 0.22 mg/ml Vs. 0.55 mg/ml). The data from Table 7 also indicates that the lower final concentration pre-lyophilization correlates with smaller change in particle size (relative to same final formulation at higher concentration). Excipient concentration ranges can include 20-100 mM Tris, 0-200 mM Arginine, 0.035-0.2% PVA10, 25-50-100-300 mM NaCl, 5-30% Sucrose and a pH of about 7.4; Reconstitution media: Normal Saline is preferred for reconstitution, however, (1) formulation buffer comprising of 50 mM HEPES, 50 mM NaCl, 9% Sucrose, 5% Glycerol pH 7.4 at 1×, 0.5× or at 0.25× concentration, (2) water, (3) D5W, (4) 0.15 M Arginine, (5) 1× or 0.5× Saline (6) 0.5 M NaCl or (7) 50 mM Tris 50 mM NaCl, 9% Sucrose, pH 7.4 at 1×, 0.5× or 0.25× may also be used in a manner such that pre-lyophilization and post-lyophilization mRNA concentration are in 1:1, 2:1 or 4:1 ratio.

TABLE-US-00008 mRNA/LNP characterization Formulation (pC-2) 50 mM HEPES, 50 mM HEPES, 50 mM Tris, 100 mM Arginine, 100 mM Arginine, 100 mM Arginine, 0.1% PVA10, 0.1% PVA10, 0.1% PVA10, 50 mM NaCl, 50 mM NaCl, 50 mM NaCl 9% Sucrose pH 7.4 9% Sucrose pH 7.4 15% Sucrose (Run1- Control) (Run 2-Control) pH 7.4 Pre- Post- Pre- Post- Pre- Post- Reconstitution Lyo Lyo Delta Lyo Lyo Delta Lyo Lyo Delta Normal Size (nm) 84 103 19 82 98 16 70.87 80.43 10 Saline (2:1 PDI 0.072 0.106 0.034 0.069 0.098 0.036 0.031 0.044 0.013 concentration) Encapsulation (%) 98.70% 90.40% -8.30% 98.90% 93.90% -4.10% 98.70% 92.00% -6.60% Normal Saline Size (nm) 84 103 19 82 97 15 70.87 82 11 (1:1 Conc.) PDI 0.072 0.088 0.016 0.069 0.086 0.024 0.031 0.065 0.034 Encapsulation (%) 98.70% 91.90% -6.80% 98.90% 94.10% -3.90% 98.70% 93.80% -4.90% Water (2:1 Size (nm) 84 130 46 NA NA conc.) PDI 0.072 0.119 0.047 NA NA Encapsulation (%) 98.70% 88.10% -10.60% NA NA 0.5M NaCl Size (nm) NA 82 99 17 NA (2:1 conc.) PDI NA 0.069 0.103 0.041 NA Encapsulation (%) NA 98.90% 94.80% -3.20% NA Size (nm) NA 82 99 17 NA 0.15M PDI NA 0.069 0.107 0.045 Arginine Encapsulation (%) NA 98.90% 94.60% -3.40% NA (2:1 conc.) D5W (2:1 Size (nm) 84 117 33 NA NA conc.) PDI 0.072 0.122 0.05 NA NA Encapsulation (%) 98.70% 90.60% -8.10% NA NA

TABLE-US-00009 Lyophilization - pC-1 Chamber Ramp Rate Pressure Step Temperature (° C.) Time (° C./min) (mTorr) Loading 20° C. NA NA NA Freezing 20° C. to -50° C. 70 min 1° C./min NA -50° C. 120 min NA NA -50° C. to -30° C. 200 min 0.1° C./min NA -30° C. 5 min NA NA Primary Drying -30° C. 2800 min NA 100 Secondary Drying -30° C. to 25° C. 550 min 0.1° C./min 100 25° C. 300 min NA 100 Storage 5° C. NA NA 100

TABLE-US-00010 Lyophilization - pC-2 Chamber Ramp Rate Pressure Step Temperature (° C.) Time (° C./min) (mTorr) Loading 18° C. 5 min NA NA Freezing 18° C. to 5° C. 13 min 1° C./min NA 5° C. 30 min NA NA 5° C. to -50° C. 45 min 1° C./min NA -50° C. 180 min NA NA Primary Drying -50° C. to -30° C. 200 min 0.1° C./min 100 -30° C. 1070 min NA 100 -30° C. to -25° C. 50 min 0.1° C./min 100 -25° C. 1200 min NA 100 Secondary Drying -25° C. to 25° C. 500 min 0.1° C./min 100 25° C. 600 min NA 100

[0101] The studies described herein examine optimal excipients and conditions that can be used to form cargo molecule/LNP within the scope of the invention described herein. Table 8 provides a list of chemicals and materials used to formulate and analyze the therapeutic cargo molecules.

TABLE-US-00011 TABLE 8 Chemical and materials list Item Vendor Catalog Number HEPES Sigma-Aldrich H3375-100G Schott borosilicate 2 cc, 13 mm serum glass vials West Pharmaceutical 68000386 Services WPS 4432/50 Grey, B2-40, stoppers West Pharmaceutical 19700004 Services Vivaspin 20, 30 kDa MWCO PES centrifugal Sartorius 26-9323-61 concentrators Slide-A-Lyzer, 20 kDa MWCO dialysis cassettes, Thermo Scientific 66003 0.5-3.0 mL Slide-A-Lyzer, 20 kDa MWCO dialysis cassettes, Thermo Scientific 66012 3-12 mL Polysorbate 80 (PS80) Croda SR40925 Poloxamer 188 (F-68) Spectrum P1169 Sucrose JT Baker 4074-05 Sorbitol EMPROVE 1.11597.9023 Potassium Sorbate Sigma-Aldrich 85520-50G Sodium

Benzotrate Sigma-Aldrich 18106-1KG-R Sodium Phosphate Monobasic Monohydrate Fisher S468-500 Trehalose Pfanstiehl T-10404 Mannitol JT Baker 2553-05 EDTA-Na2 JT Baker 8995-01 Dextran 40 Spectrum DE124 Poly(vinyl alcohol) Sigma-Aldrich 348406-25g L-Arginine Sigma-Aldrich A8094-1KG L-Glutamic Acid JT Baker M746-07 L-Histidine JT Baker 2080-06 L-Methionine Sigma-Aldrich M5308-100G Cis-4-Hydroxy-D-proline Sigma-Aldrich H5877 L-Proline Sigma-Aldrich P5607-100G L-Threonine VWR E808 Glacial Acetic Acid EMD AX0073-9 Sodium Hydroxide Macron 7680-06 BSA Sigma-Aldrich A2153 60% Iodixanol, Optiprep Density Gradient Medium Sigma-Aldrich D1556-250 mL RNase Free Water Ambion AM9932 10% TX-100 Sigma-Aldrich 93443-100 ML RiboGreen Reagent Invitrogen R11490A Ribo Ruler High Range RNA Ladder Thermo Fisher SM1821 Quant-iT™ RiboGreen™ RNA Assay Kit Invitrogen R11490 RNase Zap Thermo Fisher AM9780 3000 Series Nanosphere™ Size Standards, 100 nm Thermo Scientific 3100A

[0102] Excipients and excipient formulations were evaluated to determine how to best (1) lyophilize the LNP composition and (2) maintain a high level of fidelity between the pre-lyophilized and post-lyophilized LNP composition.

[0103] The evaluation study included various excipients to stabilize cargo molecule through lyophilization. The physical stability of cargo molecule was evaluated in three rounds of screening by analytical methods previously verified as effective stability-indicating assays for this product. These studies analyze the following: visual appearance, dynamic light scattering (DLS), and RiboGreen RNA quantitation.

[0104] For lyophilization, formulations containing 9% sucrose, 50 mM NaCl, and 50 mM Tris-HCl consistently displayed elegant, lyophilized cakes across all three rounds. RNase-free water reconstitution of sucrose formulations containing arginine and PVA10 displayed the least opalescence. Across the three rounds, DLS results following lyophilization revealed that formulations containing 9% sucrose showed smaller increases in hydrodynamic diameters compared to other non-sucrose formulations. The inclusion of arginine and PVA10 further reduced the increases in hydrodynamic diameters observed in DLS results. RiboGreen analysis concluded that reconstitution with RNase-free water yielded greater RNA percent encapsulation than reconstitution with TNS (Tris, NaCl, and sucrose) buffer. Despite experiencing smaller losses of encapsulation, iodixanol-containing formulations exhibited large increases in hydrodynamic diameters. Furthermore, despite displaying the most elegant cakes, Dextran 40 formulations were the most opalescent of the reconstituted formulations in Round 3, displayed the largest loss of encapsulation, and exhibited the greatest increase in hydrodynamic diameters. Notably, formulations containing arginine, PVA10, and NaCl showed high percentages of RNA encapsulation after reconstitution, reinforcing their beneficial effects on stability as observed in DLS and visual appearance.

[0105] In summary, the formulation of 9% sucrose, 50 mM NaCl, 50 mM Tris-HCl, 100 mM arginine, and 0.05% PVA10 presents superior stabilization of test sample through lyophilization and subsequent reconstitution with water, compared to other formulations tested. The stabilization benefits of these excipients may be further corroborated with the analysis of test sample's biological activity.

[0106] An Excipient Screening Study for test sample was conducted to assess excipients that stabilize test sample at 0.25 mg/mL during lyophilization. The active pharmaceutical ingredient (API) examined in this study was an mRNA cargo molecule expressing a GOI. The material used for this study was comprised of the following: [0107] Excipient Screening, Round 1, Drug Product (DP), 1 mg/mL in 50 mM HEPES, 50 mM NaCl, 9%, Sucrose, 5% Glycerol, pH 7.4, Lot 1119-A. [0108] Excipient Screening, Round 2 Drug Product (DP), Stored at -70° C., 0.99 mg/mL in 50 mM HEPES, 50 mM NaCl, 9% Sucrose, 5% Glycerol, pH 7.4, Lot 210113. [0109] Excipient Screening, Round 3 Drug Product (DP), Stored at -70° C., 1.00 mg/mL in 50 mM HEPES, 50 mM NaCl, 9% Sucrose, 5% Glycerol, pH 7.4, Lot 210825, Stored at -70° C.

[0110] In the Excipient screening study, the following parameters were fixed: [0111] (1) Fill volume: 0.1 mL/well in a 96-well plate [0112] (2) API Concentration: 0.25 mg/mL [0113] (3) pH: 7.4

[0114] The following excipients were examined in each round of the screening (Tables 9-11):

[0115] Round 1: HEPES, Tris-HCl, polymers, hydrophobic salts, amino acids, NaCl, sucrose, and iodixanol. [0116] Round 2: Various concentrations of PVA10, P188, arginine, proline, sucrose, and iodixanol. [0117] Round 3: EDTA, methionine, Dextran 40, various concentrations of Tris-HCl, arginine, PVA10, sucrose, and NaCl.

TABLE-US-00012 TABLE 9 Excipient Screening Matrix (Round 1) Buffer Hydrophobic salts/AA Form # (50 mM) Polymers (100 mM) NaCl Sucrose/Iodixanol 1 HEPES 0.05% Sodium Benzoate 50 mM None 2 PVA10 Potassium Sorbate 3 0.2% Sodium Benzoate 4 P188 Potassium Sorbate 5 Tris-HCl 0.05% Potassium Sorbate 6 PVA10 Sodium Benzoate 7 0.2% Potassium Sorbate 8 P188 Sodium Benzoate 9 HEPES 0.05% Potassium Sorbate None PVA10 10 0.2% Sodium Benzoate P188 11 Tris-HCl 0.05% Sodium Benzoate PVA10 12 0.2% Potassium Sorbate 13 P188 None 50 mM 2% Sucrose 14 9% Sucrose 15 0.05% PVA10 16 0.2% Arginine None 17 P188 100 mM Arginine + 100 mM Potassium Sorbate 18 100 mM Proline + 100 mM Sodium Benzoate 19 HEPES 0.05% Proline None 9% Iodixanol 20 PVA10 Arginine 21 0.2% Proline 9% Sucrose 22 P188 Arginine 23 Tris-HCl 0.05% Proline 9% Iodixanol 24 PVA10 Arginine 25 0.2% Proline 9% Sucrose 26 P188 Arginine 27 Proline 50 mM 2% Sucrose 28 20% Sucrose 29 Arginine 2% Sucrose 30 9% Sucrose

TABLE-US-00013 TABLE 10 Excipient Screening Matrix (Round 2) Buffer PVA10 P188 Amino Amino Acid Sucrose/ NaCl Form # (50 mM) (%) (%) Acid Conc. (M) Iodixanol (9%) (mM) 1 Tris-HCl 0.1 0.4 Arginine 0.05 Sucrose NA 2 Tris-HCl 0.1 0.2 Proline 0.05 Sucrose NA 3 Tris-HCl 0.05 0.4 Proline 0.1 Sucrose NA 4 Tris-HCl 0.05 NA Arginine 0.1 Sucrose NA 5 Tris-HCl 0.035 0.2 Proline 0.1 Sucrose NA 6 Tris-HCl 0.035 NA Proline 0.1 Sucrose NA 7 Tris-HCl 0.1 0.2 Arginine 0.05 Iodixanol NA 8 Tris-HCl 0.1 NA Proline 0.1 Iodixanol NA 9 Tris-HCl 0.05 0.4 Arginine 0.1 Iodixanol NA 10 Tris-HCl 0.05 0.2 Proline 0.05 Iodixanol NA 11 Tris-HCl 0.035 0.4 Proline 0.05 Iodixanol NA 12 Tris-HCl 0.035 NA Arginine 0.05 Iodixanol NA 13 Tris-HCl 0.1 0.4 Arginine 0.05 Sucrose 50 14 Tris-HCl 0.1 0.2 Proline 0.05 Sucrose 50 15 Tris-HCl 0.05 0.4 Proline 0.1 Sucrose 50 16 Tris-HCl 0.05 NA Arginine 0.1 Sucrose 50 17 Tris-HCl 0.035 0.2 Proline 0.1 Sucrose 50 18 Tris-HCl 0.035 NA Proline 0.1 Sucrose 50 19 Tris-HCl 0.1 0.2 Arginine 0.05 Iodixanol 50 20 Tris-HCl 0.1 NA Proline 0.1 Iodixanol 50 21 Tris-HCl 0.05 0.4 Arginine 0.1 Iodixanol 50 22 Tris-HCl 0.05 0.2 Proline 0.05 Iodixanol 50 23 Tris-HCl 0.035 0.4 Proline 0.05 Iodixanol 50 24 Tris-HCl 0.035 NA Arginine 0.05 Iodixanol 50 25 Tris-HCl NA 0.2 Proline 0.1 Iodixanol 50 26 Tris-HCl NA 0.2 Arginine 0.1 Iodixanol 50 27 Tris-HCl NA NA Proline 0.1 Iodixanol 50 28 Tris-HCl NA NA Arginine 0.1 Sucrose 50 29 Tris-HCl NA NA NA NA Iodixanol 50 30 Tris-HCl NA NA NA NA Sucrose 50

TABLE-US-00014 TABLE 11 Excipient Screening Matrix (Round 3) Tris- HCl Arginine PVA10 NaCl EDTA Methionine Dextran Form # (mM) (M) (%) Sucrose (%) (mM) (mM) (mM) 40 (%) 1 50 0 0 9 50 0 0 0 2 50 0 0 3 50 0 0 0 3 100 0 0 9 50 0 0 0 4 200 0 0 9 50 0 0 0 5 50 0 0 9 150 0 0 0 6 50 0 0 9 300 0 0 0 7 50 0.1 0 9 50 0 0 0 8 50 0.2 0 9 50 0 0 0 9 50 0.1 0.05 9 50 0 0 0 10 50 0.2 0.05 9 50 0 0 0 11 200 0.1 0.05 9 50 0 0 0 12 100 0.1 0.05 9 50 0 0 0 13 50 0.1 0.05 9 150 0 0 0 14 50 0.01 0.05 9 300 0 0 0 15 50 0.01 0.05 3 50 0 0 0 16 50 0 0 9 50 0.1 50 1 17 50 0.1 0.05 9 50 0.1 50 1 18 50 0 0 9 50 0.1 50 0 19 50 0.1 0.05 9 50 0.1 50 0 20 50 0 0 9 50 0.1 0 0 21 50 0.1 0.05 9 50 0.1 0 0 22 Positive Control Buffer 23 50 0 0 9 50 0 0 0 24 50 0 0 9 50 0 0 0 25 50 0.1 0.05 9 50 0 0 0 26 50 0.1 0.05 9 50 0 0 0 27 50 0 0 9 150 0 0 0 28 50 0 0 9 150 0 0 0 29 50 0.1 0.05 9 150 0 0 0 30 50 0.1 0.05 9 150 0 0 0

Formulation Preparation

[0118] Two vials of the drug product (DP) in Round 1 and Round 2 and one vial of DP in Round 3 were removed from -70° C. storage and thawed at ambient temperature. The DP was pooled when

across all three rounds.

Analytical Methods

[0125] Both pre-lyophilization and lyophilized samples were analyzed by the following analytical methods shown in Table 15.

[0126] Visual Inspection: Visual inspection was performed under a white light source (13W fluorescent tube) against black and white backgrounds. Digital photographs were acquired of all formulations at pre-lyophilization samples, lyophilized samples, and reconstituted samples.

TABLE-US-00018 TABLE 15 Analytical Methods Excipient Screening Study Round 1 Round 2 Round 3 Analytical Lyophilized & Lyophilized & Lyophilized & Methods Pre-Lyo Reconstituted Pre-Lyo Reconstituted Pre-Lyo Reconstituted Visual X X X X X X DLS X X X X X X RiboGreen X X X X X X

[0127] Visual Inspection: Visual inspection was performed under a white light source (13W fluorescent tube) against black and white backgrounds. Digital photographs were acquired of all formulations at pre-lyophilization samples, lyophilized samples, and reconstituted samples.

[0128] Dynamic Light Scattering (DLS): Changes in light intensity (scattered as particles diffuse by Brownian motion) were measured, and information about particle population, size, and size distribution were collected using the DLS DynaPro Plate Reader II. The standard was prepared by adding a drop of 3000 Series Nanosphere™ Size Standards (100 nm, Thermo Scientific, Cat #3100A) to a 1 mL aliquot of 1×PBS. For all sample sets, 100 nm Nanosphere size standard results met system suitability criteria of cumulants mean particle diameter (d. nm) of 100 nm±3 nm and the average polydispersity index (PDI)<0.1, as shown in Table A1. Samples were diluted to 0.1 mg/mL by diluting 10 µL of 0.25 mg/mL formulation in 240 µL of 1×PBS (pH 7.4). The DP control samples were diluted to 0.01 mg/mL by diluting 10 µL of 1 mg/mL DP control in 1 mL of 1×PBS. Both Nanosphere size standard and diluted samples were then degassed at 70 torr for 30 minutes to reduce the presence of air bubbles. After degassing, 25 µL/well of each sample was loaded into a 384-well plate in triplicate per sample. The plate was then centrifuged for 120 seconds at 2,000 rpm using a 5810 Eppendorf centrifuge. The plate reader was set to 25° C. and 10 acquisitions per sample at a rate of 5 seconds per acquisition, setting the total run time of 50 seconds/well. Particle sizes were analyzed in 1≤x≤10000 nm radii range. Particle radius distribution was categorized into Peak 1 (1-10 nm), Peak 2 (10-100 nm), Peak 3 (100-1000 nm), and Peak 4 (1000-5000 nm). All reported data is a cumulative average of these 10 acquisitions.

RiboGreen Analysis (RNA Quantification)

[0129] RNA standard (100 µg/mL in Tris-EDTA (TE) buffer from Quant-iT™ RiboGreen™ RNA Assay Kit) was diluted to 10 µg/mL with 1×TE buffer. The 10 µg/mL RNA standard stock was then diluted to 5, 4, 3.2, 2.4, 1.6, 0.8, and 0.4 µg/mL stock solutions with 1×TE buffer. The standard stock solutions were then diluted 1:1 with either 1×TE or 2% Triton X-100. During each round samples were prepared by diluting to ~50 µg/mL with 1×TE. Then, for Free RNA samples, an additional 1:10 dilution with TE buffer was performed to achieve 5 µg/mL; samples were not vortexed. For Total RNA samples, a 1:50 dilution was performed with 1% Triton X-100 to achieve 1 µg/mL. Each Total RNA sample was then vortexed for approximately 10 seconds for complete disruption. A RiboGreen working solution was prepared via 1:200 dilution of RiboGreen stock with 1×TE buffer. Following RNA standard and DP sample preparation, 100 µL of each sample was mixed with 100 µL of RiboGreen working solution and incubated in darkness at ambient temperature for 10 minutes. Two replicates were performed for RNA standards and three replicates were performed for DP samples. All plates were read on a fluorescence plate reader with excitation at 485 nm and emission at 530 nm. For all plates, standard curve R² values met the acceptance criteria of >0.975, and % RSD of the two standard replicates met the acceptance criteria of ≤20%. See for example, Tables A19-A26 below.

[0130] The lyophilization cycle profiles for each round of screening are shown in FIG. 6 through FIG. 8. For Round 1, primary drying at -20° C. took ~12 hrs to reach a differential of <3 mT of

Pirani Gauge/Capacitance Manometer (PVG/CM) indicating the completion of the primary drying process (FIG. 6). For Round 2 and Round 3, primary drying at -30°C . took 15-16 hours to reach a differential of $<3\text{ mT}$ of PVG/CM (FIG. 7 and FIG. 8). The secondary drying was completed at 20°C . for Round 1 and at 5°C . for Rounds 1 and 3.

Visual Observation

[0131] All formulations of cargo molecule were visually inspected before lyophilization (pre-lyo), after lyophilization (post-lyo), and after reconstitution for clarity, coloration, and visible precipitation in liquid samples, or for cake appearance in lyophilized samples. Photographs of the formulations in each plate were taken (FIG. 9 through FIG. 18). In all rounds of screening, pre-lyophilization samples (pre-lyo) at 0.25 mg/mL showed comparable visual appearances to each other, with slightly less opalescence than the DP (5°C .) control at 1 mg/mL .

[0132] After lyophilization, in Round 1, formulations containing 9% sucrose or iodixanol displayed elegant cakes with no or few minor bubbles. All the other formulations showed signs of collapse or displayed many bubbles. In Round 2, all 30 formulations contained 9% sucrose or iodixanol and displayed elegant cakes with no or few minor bubbles. In Round 3, formulations containing 9% sucrose, 50 mM NaCl , and 50 mM Tris-HCl displayed elegant cakes with a few minor bubbles. Increasing concentrations of NaCl (150 mM and 300 mM) or Tris-HCl (100 mM and 200 mM) resulted in highly porous or collapsed cakes. The Positive Control Buffer formulation displayed a porous lyophilized cake that was inflated from the well.

[0133] After reconstitution with water or 12.5 mM , Tris buffer containing 12.5 mM NaCl and 2.25% sucrose, $\text{pH } 7.5$ (TNS buffer), the various visual turbidities of the samples did not correlate to the lyophilized cake formations. In Round 1, the reconstituted formulations containing hydrophobic salts presented as white turbid liquids with undissolved precipitates, except for F17. Round 1 F17 formulation contained arginine with potassium sorbate and the reconstituted sample displayed a comparable opalescence level to its pre-lyophilization sample. Formulations containing 9% sugar (sucrose or iodixanol) and amino acids (arginine or proline) showed comparable opalescence levels to their pre-lyophilization samples. All other formulations that contained lower or higher sugar concentrations than 9% appeared to be more opalescent and turbid than their pre-lyophilization controls. In Round 2, samples were reconstituted with either filtered MQ water or TNS buffer. The reconstituted formulations containing 9% sucrose, with PVA10 or without P188, F1-F5, and F13-F17, showed comparable opalescence to their pre-lyophilization samples. All other formulations appeared to be more opalescent than their pre-lyophilization controls. The formulations were also less opalescent after reconstitution with TNS buffer compared to samples reconstituted with MQ water. In Round 3, all formulations reconstituted with RNase-free water were dissolved easily but increased in opalescence compared to their pre-lyophilization controls. Again, turbidity of the reconstituted formulations did not correlate with lyo-cake formation. For example, formulations containing Dextran 40, F16 and F17, were the most opalescent, but their lyophilized formulations displayed ideal cake formation.

Dynamic Light Scattering (DLS) Analysis

[0134] Pre-lyophilization and reconstituted post-lyophilization of test sample were analyzed by dynamic light scattering (DLS) to determine particle size distribution. Upon meeting the system suitability (Table A1), DLS results for each round of screening are shown in Tables A2-A8 and Tables 16 and 17.

TABLE-US-00019 TABLE A1 Dynamic Light Scattering, System Suitability Criteria Round 1 Round 2 Round 3 Measured Measured Measured Values Passing Values Passing Values Passing Criteria ($n = 3$) Values ($n = 3$) Values ($n = 3$) Values d. nm $100 \pm 3\text{ nm}$ $100 \pm 3\text{ nm}$ $100 \pm 3\text{ nm}$ ($100\text{ nm } 97.5\text{-} 99.1\text{-} 101.6\text{-}$ Standard)* 99.7 nm 100.9 nm 102.5 nm PDI 0.03 Average 0.04 Average 0.03 Average measurement measurement measurement must be less must be less must be less than 0.1 than 0.1 than 0.1

TABLE-US-00020 TABLE A2 Round 1, Pre-Lyo, DLS Average Main Peak d .Math. nm % Form

No. Form Code (n = 3) PDI d .Math. nm Intensity 100 nm 98.7 0.04 103.1 99.9-100.0 DP Control 69.8 0.14 76.6 99.4-99.9 F1 H-VBN 72.9 0.15 81.2 99.8-99.9 F2 H-VKN 76.2 0.13 85.7 99.6-99.8 F3 H-FBN 73.0 0.15 81.6 96.2-99.8 F4 H-FKN 75.5 0.14 83.9 99.1-99.9 F5 T-VKN 76.1 0.19 88.5 95.7-100.0 F6 T-VBN 77.4 0.19 85.9 99.7-99.9 F7 T-FKN 75.2 0.16 83.8 97.8-99.8 F8 T-FBN 75.3 0.20 85.3 99.3-99.9 F9 H-VK 76.2 0.14 88.7 99.9-99.9 F10 H-FB 74.8 0.18 84.3 95.3-99.8 F11 T-VB 78.3 0.22 93.6 99.0-100.0 F12 T-FK 74.6 0.17 86.2 96.5-99.2 F13 T-F2SN 74.5 0.17 82.2 97.3-99.9 F14 T-F9SN 76.0 0.16 81.9 99.8-100.0 F15 T-V9SN 74.9 0.19 85.0 99.8-99.9 F16 T-FRN 74.9 0.16 83.4 99.4-99.9 F17 T-FKRN 74.3 0.16 85.3 99.0-99.9 F18 T-FBPN 75.4 0.18 87.1 99.6-99.9 F19 H-VPI 71.7 0.15 81.0 98.0-99.9 F20 H-VRI 72.0 0.18 82.0 99.9-99.9 F21 H-FP9S 72.5 0.16 78.8 95.2-99.8 F22 H-FR9S 72.1 0.17 76.1 99.0-100.0 F23 T-VPI 74.4 0.17 83.8 99.6-99.9 F24 T-VRI 74.4 0.17 81.0 98.7-99.8 F25 T-FP9S 76.2 0.16 84.8 96.3-99.9 F26 T-FR9S 74.8 0.16 81.2 98.4-99.9 F27 T-FP2SN 74.1 0.17 83.6 92.4-99.9 F28 T-FP20SN 76.1 0.16 80.5 94.9-99.8 F29 T-FR2SN 74.8 0.18 86.4 98.4-99.9 F30 T-FR9SN 75.5 0.19 89.7 99.9-99.9

TABLE-US-00021 TABLE A3 Round 1, Post-Lyo, DLS Average Main Peak d .Math. nm % Form No. Form Code (n = 3) PDI d .Math. nm Intensity 100 nm 98.8 0.03 101.5 99.8-99.9 F1 H-VBN 188.6 0.45 182.5 94.3-99.3 F2 H-VKN 215.5 0.17 218.0 97.2-99.0 F3 H-FBN 171.7 0.35 175.1 92.0-95.8 F4 H-FKN 244.2 0.39 276.8 95.6-99.4 F5 T-VKN 146.1 0.44 138.4 82.4-91.0 F6 T-VBN 228.9 0.41 244.4 87.6-100.0 F7 T-FKN 155.6 Multi 166.1 88.7-99.8 F8 T-FBN 213.4 0.48 217.5 90.9-94.0 F9 H-VK 218.4 0.21 229.5 98.4-99.8 F10 H-FB 179.8 0.31 173.1 96.9-99.9 F11 T-VB 268.8 Multi 286.3 91.0-97.5 F12 T-FK 173.5 Multi 191.1 98.6-99.9 F13 T-F2SN 199.5 0.16 210.0 95.7-99.4 F14 T-F9SN 96.2 0.2 118.5 99.8-100.0 F15 T-V9SN 103.6 0.19 118.2 94.9-99.0 F16 T-FRN 151.8 0.33 167.9 87.1-97.2 F17 T-FKRN 112.1 0.39 134.7 83.3-99.4 F18 T-FBPN 208.7 0.24 222.5 96.4-99.4 F19 H-VPI 115.3 0.26 114.1 94.3-100.0 F20 H-VRI 122.2 Multi 140.9 55.9-97.4 F21 H-FP9S 131.3 0.2 132.7 97.3-100.0 F22 H-FR9S 96.0 0.32 99.6 96.7-99.9 F23 T-VPI 112.6 Multi 142.1 60.5-99.9 F24 T-VRI 128.8 Multi 168.3 93.6-100.0 F25 T-FP9S 111.0 0.24 137.6 96.4-100.0 F26 T-FR9S 95.5 0.26 106.7 94.0-99.9 F27 T-FP2SN 223.2 0.38 237.9 88.2-99.4 F28 T-FP20SN 200.6 0.21 235.1 99.1-100.0 F29 T-FR2SN 132.4 0.35 133.0 93.8-99.9 F30 T-FR9SN 91.9 0.27 105.4 93.8-99.6

TABLE-US-00022 TABLE A4 Round 2, Pre-Lyo, DLS Average Main Peak d .Math. nm % Sample (n = 3) PDI d .Math. nm Intensity 100 nm Standard 99.5 0.03 103.2 99.9-100.0 DP (5° C.) Control 85.4 0.16 84.0 99.4-99.6 F1 87.3 0.15 101.8 97.4-99.8 F2 87.6 0.16 101.4 98.4-98.5 F3 88.4 0.14 101.2 99.6-99.8 F4 87.9 0.15 103.8 99.9-100.0 F5 87.7 0.15 101.2 98.1-99.4 F6 88.1 0.15 101.5 98.5-99.9 F7 86.5 0.16 101.1 99.1-99.7 F8 87.5 0.16 100.4 99.5-99.7 F9 87.5 0.18 104.9 94.3-99.9 F10 87.1 0.16 104.1 99.4-99.9 F11 87.6 0.17 102.6 94.4-98.1 F12 88.1 0.17 102.5 99.0-99.6 F13 87.3 0.18 104.2 98.0-100.0 F14 86.8 0.15 101.2 98.1-99.8 F15 86.3 0.17 104.5 99.8-99.9 F16 86.8 0.15 96.6 97.5-99.6 F17 86.5 0.16 101.6 97.3-99.4 F18 85.9 0.2 103.2 93.9-99.9 F19 85.7 0.18 103.5 95.6-99.9 F20 85.7 0.19 103.6 98.9-99.8 F21 86.7 0.16 96.5 95.1-99.4 F22 86.9 0.17 103.7 99.4-100.0 F23 86.3 0.17 101.9 98.1-99.9 F24 86.4 0.16 101.3 98.2-99.9 F25 85.5 0.18 99.2 94.9-99.9 F26 85.6 0.17 100.2 98.2-98.6 F27 86 0.16 101.6 99.8-99.9 F28 86.4 0.17 103.3 99.6-100.0 F29 86.2 0.15 100.8 94.6-99.7 F30 87.3 0.16 99.4 93.8-99.8

TABLE-US-00023 TABLE A5 Round 2, Reconstituted with Milli-Q water, DLS Average Main Peak d .Math. nm % Sample (n = 3) PDI d .Math. nm Intensity 100 nm Standard 100.2 0.03 103.8 99.8-100.0 F1 93.0 0.18 113.4 97.0-99.9 F2 93.4 0.2 110.6 95.6-99.6 F3 92.4 0.17 110.0 95.5-99.9 F4 97.0 0.18 113.4 98.8-99.8 F5 93.4 0.18 110.7 99.7-99.9 F6 111.9 0.16 129.4 99.3-99.9 F7 125.2 0.49 156.2 99.9-99.9 F8 115.3 0.38 131.9 91.9-99.9 F9 124.7 0.49 161.3 97.3-99.9 F10 124.5 0.41 149.5 99.7-99.9 F11 125.5 0.47 154.0 96.2-100.0 F12 136.0 0.42 152.9 97.2-100.0 F13 90.9 0.19 111.5 99.6-99.9 F14 93.7 0.21 117.4 99.1-100.0 F15 92.2 0.18 110.6 99.4-99.9 F16 94.4 0.2 112.1 95.6-99.6 F17 92.3 0.18 110.8 98.8-99.9 F18 103.7 0.18 118.5 93.8-99.7 F19 110.9 0.38 132.0 91.9-100.0 F20 106.6 0.28 123.2 90.8-98.5 F21 119.2 0.49 166.5 99.4-99.9 F22 111.9 0.39 140.4 95.3-99.9 F23 114.5 0.42 134.3 97.8-100.0 F24 126.5 0.4 151.2 99.7-99.9 F25 120.7 0.48 151.2

99.8-100.0 F26 136.3 0.48 167.8 94.1-99.7 F27 127.8 0.43 142.0 90.4-100.0 F28 116.8 0.34 133.4
 92.7-99.9 F29 126.4 0.45 134.6 93.6-99.4 F30 110.1 0.2 127.1 98.7-100.0
 TABLE-US-00024 TABLE A6 Round 2, Reconstituted with TNS buffer, DLS Average Main Peak d .Math. nm % Sample (n = 3) PDI d .Math. nm Intensity 100 nm Standard 99.5 0.04 104.3 99.9-
 100.0 F1 91.7 0.18 100.8 97.3-99.7 F2 92.7 0.18 92.6 97.9-99.5 F3 90.1 0.18 84.4 99.5-99.8 F4
 99.3 0.18 89.1 99.3-100.0 F5 94.0 0.18 89.5 99.6-100.0 F6 109.8 0.16 115.5 98.2-100.0 F7
 112.7 0.38 123.1 97.6-99.9 F8 112.7 0.25 96.5 98.5-100.0 F9 123.9 0.41 133.8 90.0-99.9 F10
 118.5 0.34 122.1 96.8-100.0 F11 117.9 0.41 130.2 91.9-99.3 F12 134.2 0.34 109.6 98.4-99.8 F13
 90.0 0.19 90.3 99.5-99.9 F14 92.4 0.18 96.3 99.2-100.0 F15 91.8 0.17 101.2 99.0-100.0 F16 95.9
 0.2 107.8 98.8-99.6 F17 94.4 0.18 102.6 94.1-99.8 F18 106.0 0.18 114.4 96.4-98.9 F19 101.4 0.3
 106.4 90.8-99.8 F20 104.2 0.28 110.6 98.8-100.0 F21 109.2 0.43 121.4 99.6-99.9 F22 104.3 0.28
 102.4 90.5-99.0 F23 108.1 0.38 127.1 98.0-99.9 F24 122.7 0.39 142.0 99.5-99.9 F25 111.8 0.38
 121.3 95.0-99.9 F26 120.1 0.39 122.4 93.5-99.9 F27 120.5 0.44 131.0 92.5-99.9 F28 113.3 0.37
 117.1 91.4-99.9 F29 127.5 0.29 92.2 100.0-100.0 F30 105.0 0.19 116.8 95.5-99.9

TABLE-US-00025 TABLE A7 Round 3, Pre-Lyo, DLS Average Main Peak d .Math. nm % Sample
 (n = 3) PDI d .Math. nm Intensity 100 nm Standard 100.1 0.03 102.1 99.7-100.0 DP (5° C.)
 Control 87.4 0.2 96.1 95.9-99.6 FB1 89.1 0.21 100.0 99.3-99.8 FB2 85.9 0.18 90.5 99.6-99.9 FB3
 86.2 0.19 102.7 99.1-99.9 FB4 89.9 0.25 97.3 99.9-99.9 FB5 86.0 0.19 98.9 99.7-99.9 FB6 85.5
 0.17 98.4 98.9-99.9 FB7 86.2 0.18 97.3 97.1-99.9 FB8 89.2 0.22 94.2 97.6-99.8 FB9 86.2 0.2 100.6
 99.7-99.9 FB10 90.3 0.2 98.4 97.5-99.9 FB11 88.1 0.2 101.8 99.3-100.0 FB12 86.3 0.19 102.0
 99.5-99.9 FB13 87.7 0.2 99.8 99.7-99.9 FB14 85.9 0.19 100.0 98.5-99.9 FB15 88.1 0.19 101.7
 98.3-99.9 FB16 86.7 0.2 97.1 99.6-99.9 FB17 88.3 0.19 108.6 99.8-100.0 FB18 85.8 0.18 94.2
 94.0-99.3 FB19 87.0 0.22 103.1 99.5-99.9 FB20 85.3 0.17 94.5 99.7-99.9 FB21 84.8 0.19 99.4
 99.6-99.9 FB22 82.7 0.19 94.9 98.1-99.6 FB23 85.2 0.18 96.4 99.8-99.9 FB24 84.2 0.19 89.8 99.0-
 99.9 FB25 85.3 0.18 98.4 98.1-100.0 FB26 85.0 0.18 101.9 99.8-100.0 FB27 87.9 0.18 103.8
 99.5-99.8 FB28 84.5 0.16 95.0 92.9-99.8 FB29 85.9 0.18 93.7 98.8-99.9 FB30 84.9 0.19 100.1
 98.6-100.0

TABLE-US-00026 TABLE A8 Round 3, Reconstituted with RNase-free water, DLS Average Main
 Peak d .Math. nm % Sample (n = 3) PDI d .Math. nm Intensity 100 nm Standard 100.1 0.03 .1
 99.7-100.0 F1 108.8 0.26 129.3 98.4-99.9 F2 176.1 0.19 217.3 90.8-99.6 F3 116.9 0.24 139.8
 99.6-100.0 F4 139.1 0.14 158.5 98.2-100.0 F5 123.8 0.29 135.4 97.5-99.9 F6 134.0 0.29 139.8
 94.1-99.9 F7 118.8 0.30 131.6 98.8-99.9 F8 145.3 0.44 160.2 80.4-99.9 F9 95.2 0.22 116.6 99.9-
 100.0 F10 105.8 0.22 115.1 95.4-99.8 F11 116.8 0.13 128.8 97.3-99.5 F12 103.2 0.22 124.9 99.6-
 99.9 F13 103.5 0.21 113.1 99.1-99.9 F14 113.4 0.17 124.4 89.9-98.3 F15 134.4 0.16 155.7 98.4-
 99.9 F16 208.5 0.29 203.6 97.2-99.6 F17 249.0 0.46 271.3 94.4-99.5 F18 111.6 0.26 112.1 98.6-
 99.9 F19 105.0 0.27 111.3 93.5-99.9 F20 119.8 0.25 141.2 99.5-99.9 F21 97.7 0.24 120.2 99.8-99.9
 F22 189.5 0.22 223.8 97.4-99.7 F23 127.3 0.27 161.1 99.4-99.9 F24 121.7 0.24 154.8 99.7-99.9
 F25 95.2 0.23 108.8 98.5-99.8 F26 93.9 0.19 108.1 88.2-99.5 F27 127.2 0.25 155.7 99.0-99.9 F28
 124.9 0.26 152.0 98.1-100.0 F29 111.5 0.14 123.9 99.7-99.9 F30 99.3 0.23 111.4 96.9-100.0

[0135] All pre-lyophilization formulations showed comparable hydrodynamic diameters and main peak intensity percentages to DP control (Table 16). Following lyophilization, however, the reconstituted formulations showed increases in hydrodynamic diameters and polydispersity (PDI) compared to their pre-lyophilization controls (Table 17).

[0136] In Round 1, formulations containing 9% sucrose, formulations containing PVA10 with sucrose or iodixanol, and F17 showed smaller increases in hydrodynamic diameters. Among them, formulations containing 9% sucrose with either poloxamer 188, arginine, or salt showed smaller increases in hydrodynamic diameters compared to the other formulations. In Round 2, reconstitution with MQ or TNS (pH 7.5) following lyophilization yielded similar results. All formulations showed increases in hydrodynamic diameters compared to their pre-lyophilization controls. Formulations containing 9% sucrose, showed smaller increases in hydrodynamic

diameters compared to those containing iodixanol. Formulations without P188 showed larger increases than their P188-containing counterparts. All formulations showed comparable main peak intensity percentages. In Round 3, all formulations were reconstituted with RNase-free water. All reconstituted formulations showed increases in hydrodynamic diameters, 46% to 150%, compared to those in their pre-lyophilization state. Formulations containing Dextran 40 showed the largest increases in hydrodynamic diameters. Formulations containing arginine and PVA10 showed smaller increases in hydrodynamic diameters. All formulation samples showed comparable main peak intensity percentages.

TABLE-US-00027 TABLE 16 Dynamic Light Scattering, Pre-lyophilization formulations Average Main Peak diameter Cumulative Main Peak Round (nm) Diameter Intensity No. Formulation (n = 3) PDI (nm) (%) Round 1 DP Control 69.80 0.14 76.60 99.4-99.9 Ave. 74.80 0.17 84.08 95.2-100.0 Formulations Round 2 DP Control 85.40 0.16 84.00 99.4-99.6 Ave. 86.85 0.16 101.76 93.8-100.0 Formulations Round 3 DP Control 87.10 0.20 96.10 95.9-99.6 Ave. 86.47 0.19 98.48 92.9-100.0 Formulations

TABLE-US-00028 TABLE 17 Dynamic Light Scattering Top Performing Formulations of Each Round Cumulative Main Peak Main Average Cumulative Peak Round Form. Diameter Diameter Intensity No. No. (nm) PDI (nm) (%) Round 1 Lyo: Milli-Q filtered water Reconstituted F22 96 0.32 99.6 96.7-99.9 F30 91.9 0.27 105.4 93.8-99.6 F26 95.5 0.26 106.7 94.0-99.9 F19 115.3 0.26 114.1 94.3-100.0 F15 103.6 0.19 118.2 94.9-99.0 F14 96.2 0.2 118.5 99.8-100.0 Round 2 Lyo: Milli-Q filtered water Reconstituted F3 92.4 0.17 110 95.5-99.9 F2 93.4 0.2 110.6 95.6-99.6 F15 92.2 0.18 110.6 99.4-99.9 F5 93.4 0.18 110.7 99.7-99.9 F17 92.3 0.18 110.8 98.8-99.9 F13 90.9 0.19 111.5 99.6-99.9 F16 94.4 0.2 112.1 95.6-99.6 F1 93 0.18 113.4 97.0-99.9 F4 97 0.18 113.4 98.8-99.8 Lyo: TNS, pH 7.5, Reconstituted F1 91.7 0.18 84.4 97.3-99.7 F2 92.7 0.18 89.1 97.9-99.5 F3 90.1 0.18 89.5 99.5-99.8 F4 99.3 0.18 90.3 99.3-100.0 F5 94 0.18 92.2 99.6-100.0 F6 109.8 0.16 92.6 98.2-100.0 Round 3 Lyo: RNase-Free Water Reconstituted F26 93.9 0.19 108.1 88.2-99.5 F25 95.2 0.23 108.8 98.5-99.8 F19 105 0.27 111.3 93.5-99.9 F30 99.3 0.23 111.4 96.9-100.0 F18 111.6 0.26 112.1 98.6-99.9 F13 103.5 0.21 113.1 99.1-99.9 F10 105.8 0.22 115.1 95.4-99.8 * n = 3 ** Ranked by Main Peak Cumulative Diameter, from smallest to largest

RiboGreen Analysis

[0137] In each round of this study, pre-lyophilization and post-lyophilization formulations of test sample were analyzed with RiboGreen RNA-binding fluorescence assay to quantify RNA concentrations. mRNA concentrations of the formulations before and after treatment with TritonX-100 (TX100) were used to calculate percentages of encapsulated mRNA. All assay plates met system suitability; results can be found in Tables A17-A26. Complete RiboGreen assay results for each round of experiments are shown in Tables A8-A16.

[0138] In each round of the study, pre-lyophilized samples of all formulations showed comparable percent encapsulation to the DP control (Table 18). Lyophilization and reconstitution with MQ water in Round 1 resulted in decreased encapsulation, with losses ranging from -9.4% to -74.4%, compared to pre-lyophilization samples. Following lyophilization, PVA-containing formulations with iodixanol or sucrose showed higher encapsulation percentages compared to the other formulations, as is seen in Table 19. The addition of arginine or NaCl also resulted in some retention of encapsulation during lyophilization.

[0139] In Round 2, lyophilization and reconstitution with MQ water or TNS buffer led to decreased encapsulation compared to pre-lyophilization samples. Losses for reconstitution with MQ water ranged from -7.1% to -61.1%, and reconstitution with TNS buffer ranged from -9.9% to -70.2%. However, reconstitution with MQ water yielded better encapsulation results than reconstitution with TNS buffer for all 30 formulations (+0.3% to +12.2%). Overall, reconstitution with MQ water is preferable compared to TNS buffer. Formulations containing iodixanol experienced smaller losses of encapsulated RNA after reconstitution, for both reconstitutions with MQ water and TNS buffer. The addition of NaCl continued to show benefits in preserving encapsulation through

lyophilization. The top-performing formulations for Round 2 can be found in Table 19. [0140] In Round 3, lyophilization and reconstitution with RNase-free water resulted in decreased percent encapsulation compared to pre-lyophilization samples. Losses ranged from -3.8% to -42.7%. Formulations containing PVA10 and arginine experienced the smallest losses of encapsulated RNA after lyophilization. High concentrations of Tris-HCl and NaCl also resulted in smaller losses of encapsulated RNA after lyophilization. Both EDTA and methionine appeared to have no significant effect in either direction, while formulations with Dextran 40 had the greatest loss of encapsulated RNA after lyophilization. Overall, the groups of F9, F25, and F26 and F13, F29, and F30 containing PVA10 and arginine experienced smaller losses of encapsulated RNA, as seen in Table 19.

TABLE-US-00029 TABLE 18 Pre- Lyophilization Formulations of Each Round by RiboGreen assay

Free RNA	Total RNA	Round	Form.	Conc.	Conc.	% No.	No.	Formulation	(ug/mL)	(ug/mL)	Encap.
1 DP Control	76.8	1208.9	93.6	Ave. Formulations	16.9	241.4	93.0	2 DP Control	10.4	1145.8	99.1
Ave. Formulations	5.3	273.5	98.1	3 DP Control	11.4	980.6	98.8	Ave. Formulations	3.1	240.1	98.7

TABLE-US-00030 TABLE 19 Lyophilization Formulations of Each Round by RiboGreen assay

Free RNA	Total RNA	Round	Form.	Conc.	Conc.	% No.	No.	No.	(μg/mL)	(μg/mL)	Encapsulation
Round 1 Lyo: Milli-Q Reconstituted	F15	37.8	241.1	84.3	F24	41.2	248.2	83.4	F20	39.1	224.7
	82.6	F19	50.6	252.4	79.9	F23	54.4	266.9	79.6	F17	81.9
	283.1	71.1	Round 2 Lyo: Milli-Q Reconstituted	F30	19.2	235.0	91.8	F29	21.4	229.5	90.7
	F20	21.3	222.6	90.4	F24	24.0	230.6	89.6	F28	24.3	223.4
	89.1	F19	25.6	228.8	88.8	F8	29.6	247.2	88.0	Lyo: TNS, pH 7.5, Reconstituted	F29
	243.4	89.2	F16	27.4	217.2	87.4	F24	30.2	231.3	86.9	F20
	30.4	227.0	86.6	F28	29.2	213.4	86.3	F19	33.2	232.5	
	85.7	F30	33.1	220.3	85.0	Round 3 Lyo: RNase-Free Water Reconstituted	F14	11.3	216.9	94.8	F29
	10.3	191.3	94.6	F11	12.5	198.1	93.7	F13	15.5	214.1	92.7
	F30	15.2	201.9	92.5	F10	18.1	221.7	91.8	F4	18.7	
	193.3	90.3	F12	21.6	210.6	89.7	F26	20.9	198.2	89.5	F15
	23.9	224.5	89.4	TABLE-US-00031 TABLE A9 Round 1, Pre-Lyo, RiboGreen assay	Free RNA	Total RNA	Form	Form	Conc.	Conc.	No.

Form	Conc.	Conc.	No.	Code	(μg/mL)	(μg/mL)	% Encap.					
DP Control	76.8	1208.9	93.6	F1 H-VBN	15.1	228.6	93.4					
F2 H-VKN	17.8	239.0	92.6	F3 H-FBN	15.8	231.1	93.2					
F4 H-FKN	20.6	251.3	91.8	F5 T-VKN	18.6	247.0	92.5					
F6 T-VBN	18.1	248.2	92.7	F7 T-FKN	14.9	239.0	93.8					
F8 T-FBN	19.6	251.9	92.2	F9 H-VK	20.3	266.1	92.4					
F10 H-FB	18.9	245.0	92.3	F11 T-VB	13.5	231.0	94.2					
F12 T-FK	12.5	227.2	94.5	F13 T-F2SN	16.4	239.3	93.1					
F14 T-F9SN	17.4	232.9	92.5	F15 T-V9SN	14.9	236.1	93.7					
F16 T-FRN	14.2	234.8	93.9	F17 T-FKRN	20.7	262.0	92.1					
F18 T-FBPN	19.3	240.7	92.0	F19 H-VPI	14.0	227.0	93.8					
F20 H-VRI	12.9	223.2	94.2	F21 H-FP9S	16.9	235.1	92.8					
F22 H-FR9S	17.8	232.0	92.3	F23 T-VPI	16.0	232.6	93.1					
F24 T-VRI	14.8	232.0	93.6	F25 T-FP9S	18.8	260.1	92.8					
F26 T-FR9S	15.7	247.5	93.6	F27 T-FP2SN	18.0	252.5	92.9					
F28 T-FP20SN	17.5	247.5	92.9	F29 T-FR2SN	17.5	247.5	92.9					
F30 T-FR9SN	18.9	254.4	92.6	TABLE-US-00032 TABLE A10 Round 1, Reconstituted with Milli-Q water, RiboGreen assay	1	Free RNA	Total RNA	Form	Form	Conc.	Conc.	No.

Code	(μg/mL)	(μg/mL)	% Encap.
F1 H-VBN	191.9	248.2	22.7
F2 H-VKN	185.6	258.5	28.2
F3 H-FBN	187.0	245.8	23.9
F4 H-FKN	162.3	204.7	20.7
F5 T-VKN	164.8	209.0	21.1
F6 T-VBN	166.9	250.6	33.4
F7 T-FKN	172.5	245.8	29.8
F8 T-FBN	172.2	253.0	31.9
F9 H-VK	156.9	193.8	19.1
F10 H-FB	165.7	221.1	25.0
F11 T-VB	168.4	226.5	25.6
F12 T-FK	146.6	183.5	20.1
F13 T-F2SN	148.2	235.6	37.1
F14 T-F9SN	118.1	246.5	52.1
F15 T-V9SN	37.8	241.1	84.3
F16 T-FRN	114.9	286.5	59.9
F17 T-FKRN	81.9	283.1	71.1
F18 T-FBPN	167.6	269.3	37.8
F19 H-VPI	50.6	252.4	79.9
F20 H-VRI	39.1	224.7	82.6
F21 H-FP9S	145.7	254.2	42.7
F22 H-FR9S	119.0	253.0	53.0
F23 T-VPI	54.4	266.9	79.6
F24 T-VRI	41.2	248.2	83.4
F25 T-FP9S	140.1	281.2	50.2
F26 T-FR9S	124.3	275.2	54.8
F27 T-FP2SN	171.4	262.6	34.7
F28 T-FP20SN	123.1	217.7	43.5
F29 T-FR2SN	91.1	251.8	63.8
F30 T-FR9SN	113.7	247.6	54.1

TABLE-US-00033 TABLE A11 RiboGreen, Pre-lyophilization Samples, Round 2 Free RNA Total RNA Conc. Conc. Sample (μg/mL) (μg/mL) % Encap. DP (5° C.) 10.4 1145.8 99.1 Control F1 6.9 285.6 97.6 F2 7.2 271.9 97.3 F3 7.2 270.6 97.3 F4 7.1 275.8 97.4 F5 6.6 273.8 97.6 F6 6.8 266.7 97.4 F7 6.6 273.8 97.6 F8 6.9 277.7 97.7 F9 5.9 277.5 97.9 F10 5.4 275.5 98.0 F11 5.4 263.5 97.9

F12 4.9 281.5 98.2 F13 5.6 274.8 98.0 F14 5.6 274.2 98.0 F15 5.8 269.5 97.9 F16 6.0 274.8 97.8
F17 1.8 264.7 99.3 F18 4.9 266.1 98.2 F19 5.1 260.7 98.0 F20 4.1 262.1 98.4 F21 4.9 265.4 98.1
F22 5.1 272.1 98.1 F23 4.6 269.4 98.3 F24 4.8 281.5 98.3 F25 4.6 278.2 98.4 F26 4.4 276.8 98.4
F27 3.4 284.4 98.8 F28 4.2 263.8 98.4 F29 3.2 282.3 98.9 F30 4.3 291.3 98.5

TABLE-US-00034 TABLE A12 RiboGreen, Lyo-milli-Q Water Reconstituted Samples, Round 2
Free RNA Total RNA Conc. Conc. Sample (µg/mL) (µg/mL) % Encap. F1 127.6 232.9 45.2 F2
144.9 239.1 39.4 F3 147.1 230.5 36.2 F4 37.3 224.5 83.4 F5 140.7 244.9 42.6 F6 62.5 243.1 74.3
F7 36.9 246.7 85.0 F8 29.6 247.2 88.0 F9 52.6 242.1 78.3 F10 46.7 247.4 81.1 F11 69.0 240.3 71.3
F12 32.2 228.9 85.9 F13 114.2 231.9 50.7 F14 125.2 240.3 47.9 F15 127.2 229.5 44.6 F16 28.6
233.7 87.8 F17 123.3 241.8 49.0 F18 46.2 217.7 78.8 F19 25.6 228.8 88.8 F20 21.3 222.6 90.4 F21
35.4 232.5 84.8 F22 30.6 229.4 86.7 F23 39.8 233.1 82.9 F24 24.0 230.6 89.6 F25 68.0 241.7 71.9
F26 59.1 232.6 74.6 F27 37.8 238.1 84.1 F28 24.3 223.4 89.1 F29 21.4 229.5 90.7 F30 19.2 235.0
91.8

TABLE-US-00035 TABLE A13 RiboGreen, Lyo-TNS Buffer, pH 7.5 Reconstituted Samples,
Round 2 Free RNA Total RNA Conc. Conc. Sample (µg/mL) (µg/mL) % Encap. F1 134.4 212.5
36.7 F2 143.8 214.3 32.9 F3 149.8 210.7 28.9 F4 35.4 208.8 83.1 F5 149.5 215.0 30.4 F6 58.6
211.3 72.2 F7 34.0 218.6 84.5 F8 37.0 234.6 84.2 F9 47.4 206.8 77.1 F10 50.6 204.4 75.3 F11 70.9
219.7 67.7 F12 33.6 209.9 84.0 F13 123.2 214.8 42.7 F14 128.2 217.2 41.0 F15 129.9 215.4 39.7
F16 27.4 217.2 87.4 F17 127.5 220.2 42.1 F18 49.1 210.4 76.7 F19 33.2 232.5 85.7 F20 30.4 227.0
86.6 F21 48.3 241.1 80.0 F22 38.3 234.3 83.6 F23 41.7 228.2 81.7 F24 30.2 231.3 86.9 F25 68.6
222.8 69.2 F26 74.2 237.8 68.8 F27 50.5 227.1 77.8 F28 29.2 213.4 86.3 F29 26.2 243.4 89.2 F30
33.1 220.3 85.0

TABLE-US-00036 TABLE A14 RiboGreen, Pre-lyophilization Samples, Round 3 Free RNA Total
RNA Conc. Conc. Sample (µg/mL) (µg/mL) % Encap. DP (5° C.) 11.4 980.6 98.8 Control FB1 3.1
227.7 98.6 FB2 3.3 231.9 98.6 FB3 3.1 225.3 98.6 FB4 3.3 218.1 98.5 FB5 2.8 231.9 98.8 FB6 3.2
233.7 98.6 FB7 3.1 237.3 98.7 FB8 3.1 227.1 98.6 FB9 3.4 248.5 98.6 FB10 3.5 242.3 98.6 FB11
4.1 246.0 98.3 FB12 3.1 226.7 98.6 FB13 3.7 246.0 98.5 FB14 3.6 251.0 98.6 FB15 3.2 247.9 98.7
FB16 3.1 244.8 98.7 FB17 3.1 244.3 98.7 FB18 3.0 229.4 98.7 FB19 2.9 248.9 98.8 FB20 3.0
247.6 98.8 FB21 2.7 255.4 98.9 FB22 4.2 254.7 98.4 FB23 3.1 265.1 98.8 FB24 3.0 248.2 98.8
FB25 2.7 238.1 98.8 FB26 2.6 240.7 98.9 FB27 2.7 231.7 98.8 FB28 2.6 229.8 98.8 FB29 2.8
241.3 98.8 FB30 3.0 242.6 98.7

TABLE-US-00037 TABLE A15 RiboGreen, Lyo-RNase-Free Water Reconstituted Samples,
Round 3 Free RNA Total RNA Conc. Conc. % Encap. Sample (µg/mL) (µg/mL) % Encap. Rank
FB1 44.5 177.6 75.0 27 FB2 16.7 112.7 85.2 17 FB3 36.8 200.0 81.6 20 FB4 18.7 193.3 90.3 7
FB5 42.9 190.3 77.5 25 FB6 34.6 187.9 81.6 21 FB7 25.3 196.4 87.1 13 FB8 41.4 196.4 78.9 24
FB9 27.8 218.2 87.3 12 FB10 18.1 221.7 91.8 6 FB11 12.5 198.1 93.7 3 FB12 21.6 210.6 89.7 8
FB13 15.5 214.1 92.7 4 FB14 11.3 216.9 94.8 1 FB15 23.9 224.5 89.4 10 FB16 79.5 180.7 56.0 30
FB17 73.2 183.1 60.0 29 FB18 45.1 193.3 76.6 26 FB19 28.0 193.3 85.5 16 FB20 39.9 212.3 81.2
23 FB21 22.3 197.4 88.7 11 FB22 56.4 187.3 69.9 28 FB23 30.5 219.5 86.1 14 FB24 31.8 205.8
84.5 18 FB25 27.7 194.5 85.7 15 FB26 20.9 198.2 89.5 9 FB27 32.5 188.9 82.8 19 FB28 35.8
192.0 81.4 22 FB29 10.3 191.3 94.6 2 FB30 15.2 201.9 92.5 5

TABLE-US-00038 TABLE A16 RiboGreen, Lyo-RNase-Free Water Reconstituted Samples,
Ranked by % Encapsulation, Round 3 Free RNA Total RNA Conc. Conc. % Encap. Sample
(µg/mL) (µg/mL) % Encap. Rank FB14 11.3 216.9 94.8 1 FB29 10.3 191.3 94.6 2 FB11 12.5 198.1
93.7 3 FB13 15.5 214.1 92.7 4 FB30 15.2 201.9 92.5 5 FB10 18.1 221.7 91.8 6 FB4 18.7 193.3
90.3 7 FB12 21.6 210.6 89.7 8 FB26 20.9 198.2 89.5 9 FB15 23.9 224.5 89.4 10 FB21 22.3 197.4
88.7 11 FB9 27.8 218.2 87.3 12 FB7 25.3 196.4 87.1 13 FB23 30.5 219.5 86.1 14 FB25 27.7 194.5
85.7 15 FB19 28.0 193.3 85.5 16 FB2 16.7 112.7 85.2 17 FB24 31.8 205.8 84.5 18 FB27 32.5
188.9 82.8 19 FB3 36.8 200.0 81.6 20 FB6 34.6 187.9 81.6 21 FB28 35.8 192.0 81.4 22 FB20 39.9
212.3 81.2 23 FB8 41.4 196.4 78.9 24 FB5 42.9 190.3 77.5 25 FB18 45.1 193.3 76.6 26 FB1 44.5

177.6 75.0 27 FB22 56.0 187.3 69.9 28 FB17 73.2 183.1 60.0 29 FB16 79.5 180.7 56.0 30

TABLE-US-00039 TABLE A17 RiboGreen, Standard Curves, Pre-lyophilization Samples Round 1
 Accep- Mea- Accep- Mea- Standard tance sured Standard tance sured Set 1 Criteria Values Set 2
 Criteria Values Standard >0.975 >0.995 Standard >0.975 >0.995 Curve R.sup.2 Curve R.sup.2
 Standard % ≤20% <2.9% Standard % ≤20% <5.1% RSD of 2 RSD of 2 replicates replicates
 Standard TE Curve: y = Standard TE Curve: y = Curve 296.54x + 4.5676 Curve 301.69x + 4.0823
 Equation TX Curve: y = Equation TX Curve: y = 271.91x + 1.3341 261.25x + 3.1302 Accep- Mea-
 Accep- Mea- Standard tance sured Standard tance sured Set 3 Criteria Values Set 4 Criteria Values
 Standard >0.975 >0.990 Standard >0.975 >0.996 Curve R.sup.2 Curve R.sup.2 Standard % ≤20%
 <5.2% Standard % ≤20% <2.7% RSD of 2 RSD of 2 replicates replicates Standard TE Curve: y =
 Standard TE Curve: y = Curve 283.35x + 7.9903 Curve 296.81x + 3.484 Equation TX Curve: y =
 Equation TX Curve: y = 266.85x + 3.526 264.93x + 2.1919

TABLE-US-00040 TABLE A18 RiboGreen, Standard Curves, Lyophilization Samples Round 1
 Accep- Mea- Accep- Mea- Standard tance sured Standard tance sured Set 1 Criteria Values Set 2
 Criteria Values Standard >0.975 >0.999 Standard >0.975 >0.999 Curve R.sup.2 Curve R.sup.2
 Standard % ≤20% <2.6% Standard % ≤20% <2.7% RSD of 2 RSD of 2 replicates replicates
 Standard TE Curve: y = Standard TE Curve: y = Curve 316.44x - 1.8744 Curve 319.52x - 0.8618
 Equation TX Curve: y = Equation TX Curve: y = 276.14x - 0.9037 275.25x - 0.5433 Accep- Mea-
 Accep- Mea- Standard tance sured Standard tance sured Set 3 Criteria Values Set 4 Criteria Values
 Standard >0.975 >0.999 Standard >0.975 >0.998 Curve R.sup.2 Curve R.sup.2 Standard % ≤20%
 <2.7% Standard % ≤20% <4.5% RSD of 2 RSD of 2 replicates replicates Standard TE Curve: y =
 Standard TE Curve: y = Curve 322.68x - 2.0813 Curve 321.27x + 1.998 Equation TX Curve: y =
 Equation TX Curve: y = 276.51x - 0.9131 278.21x - 0.2791

TABLE-US-00041 TABLE A19 RiboGreen, Standard Curves, Pre-lyophilization Samples Round 2
 Accep- Mea- Accep- Mea- Standard tance sured Standard tance sured Set 1 Criteria Values Set 2
 Criteria Values Standard >0.975 >0.997 Standard >0.975 >0.996 Curve R.sup.2 Curve R.sup.2
 Standard % ≤20% <2.9% Standard % ≤20% <3.9% RSD of 2 RSD of 2 replicates replicates
 Standard TE Curve: y = Standard TE Curve: y = Curve 337.03x + 0.6133 Curve 328.13x + 3.453
 Equation TX Curve: y = Equation TX Curve: y = 255.63x + 2.999 251.11x + 5.6462 Accep- Mea-
 Accep- Mea- Standard tance sured Standard tance sured Set 3 Criteria Values Set 4 Criteria Values
 Standard >0.975 >0.996 Standard >0.975 >0.991 Curve R.sup.2 Curve R.sup.2 Standard % ≤20%
 <5.4% Standard % ≤20% <5.1% RSD of 2 RSD of 2 replicates replicates Standard TE Curve: y =
 Standard TE Curve: y = Curve 322.29x + 4.8151 Curve 321.03x + 7.3171 Equation TX Curve: y =
 Equation TX Curve: y = 249.2x + 4.5573 242.45x + 7.1056

TABLE-US-00042 TABLE A20 RiboGreen, Standard Curves, Lyo-MQ Reconstituted Samples
 Round 2 Accep- Mea- Accep- Mea- Standard tance sured Standard tance sured Set 1 Criteria
 Values Set 2 Criteria Values Standard >0.975 >0.997 Standard >0.975 >0.997 Curve R.sup.2 Curve
 R.sup.2 Standard % ≤20% <3.9% Standard % ≤20% <7.2% RSD of 2 RSD of 2 replicates
 replicates Standard TE Curve: y = Standard TE Curve: y = Curve 363.7x - 1.3867 Curve 372.41x -
 1.9977 Equation TX Curve: y = Equation TX Curve: y = 278.73x + 1.002 278.01x + 1.0819 Accep-
 Mea- Accep- Mea- Standard tance sured Standard tance sured Set 3 Criteria Values Set 4 Criteria
 Values Standard >0.975 >0.998 Standard >0.975 >0.996 Curve R.sup.2 Curve R.sup.2 Standard %
 ≤20% <5.5% Standard % ≤20% <7.6% RSD of 2 RSD of 2 replicates replicates Standard TE
 Curve: y = Standard TE Curve: y = Curve 369.48x + 0.5963 Curve 359.68x + 5.5513 Equation TX
 Curve: y = Equation TX Curve: y = 269.55x + 0.9963 273.01x + 2.6739

TABLE-US-00043 TABLE A21 RiboGreen, Standard Curves, Lyo-TNS Buffer, pH 7.5
 Reconstituted Samples Round 2 Accep- Mea- Accep- Mea- Standard tance sured Standard tance
 sured Set 1 Criteria Values Set 2 Criteria Values Standard >0.975 >0.998 Standard >0.975 >0.998
 Curve R.sup.2 Curve R.sup.2 Standard % ≤20% <3.6% Standard % ≤20% <3.1% RSD of 2 RSD of
 2 replicates replicates Standard TE Curve: y = Standard TE Curve: y = Curve 344.83x - 5.5663

Curve $343.51x - 4.7192$ Equation TX Curve: $y = 271.19x - 0.5859$
 272.74x - 1.4877 Accep- Mea- Accep- Mea- Standard tance sured Standard tance sured Set 3
 Criteria Values Set 4 Criteria Values Standard >0.975 >0.998 Standard >0.975 >0.998 Curve
 R.sup.2 Curve R.sup.2 Standard % $\leq 20\%$ <5.1% Standard % $\leq 20\%$ <5.6% RSD of 2 RSD of 2
 replicates replicates Standard TE Curve: $y =$ Standard TE Curve: $y =$ Curve $347.54x - 3.9121$
 Curve $348.21x - 1.4017$ Equation TX Curve: $y =$ Equation TX Curve: $y = 271.93x - 0.1143$
 $265.92x + 2.5283$

TABLE-US-00044 TABLE A22 RiboGreen, Standard Curves, Pre-lyophilization Samples Round 3
 Accep- Mea- Accep- Mea- Standard tance sured Standard tance sured Set 1 Criteria Values Set 2
 Criteria Values Standard >0.975 >0.997 Standard >0.975 >0.998 Curve R.sup.2 Curve R.sup.2
 Standard % $\leq 20\%$ <4.1% Standard % $\leq 20\%$ <4.7% RSD of 2 RSD of 2 replicates replicates
 Standard TE Curve: $y =$ Standard TE Curve: $y =$ Curve $342.37x - 0.9123$ Curve $346.26x + 0.8912$
 Equation TX Curve: $y =$ Equation TX Curve: $y = 278.32x + 1.5895$ $267.47x + 1.5509$ Accep- Mea-
 Accep- Mea- Standard tance sured Standard tance sured Set 3 Criteria Values Set 4 Criteria Values
 Standard >0.975 >0.998 Standard >0.975 >0.998 Curve R.sup.2 Curve R.sup.2 Standard % $\leq 20\%$
 <4.8% Standard % $\leq 20\%$ <5.8% RSD of 2 RSD of 2 replicates replicates Standard TE Curve: $y =$
 Standard TE Curve: $y =$ Curve $340.32x - 0.1105$ Curve $336.42x + 1.086$ Equation TX Curve: $y =$
 Equation TX Curve: $y = 257.16x + 2.6614$ $260.79x + 1.9737$

TABLE-US-00045 TABLE A23 RiboGreen, Standard Curves, Lyo-RNase-Free Water
 Reconstituted Samples Round 3 Accep- Mea- Accep- Mea- Standard tance sured Standard tance
 sured Set 1 Criteria Values Set 2 Criteria Values Standard >0.975 >0.994 Standard >0.975 >0.994
 Curve R.sup.2 Curve R.sup.2 Standard % $\leq 20\%$ <3.1% Standard % $\leq 20\%$ <17.7% RSD of 2 RSD
 of 2 replicates replicates Standard TE Curve: $y =$ Standard TE Curve: $y =$ Curve $351.75x + 4.9211$
 Curve $342.82x + 1.7801$ Equation TX Curve: $y =$ Equation TX Curve: $y = 274.97x + 6.6726$
 $240.04x + 9.7258$ Accep- Mea- Accep- Mea- Standard tance sured Standard tance sured Set 3
 Criteria Values Set 4 Criteria Values Standard >0.975 >0.995 Standard >0.975 >0.991 Curve
 R.sup.2 Curve R.sup.2 Standard % $\leq 20\%$ <2.7% Standard % $\leq 20\%$ <10.7% RSD of 2 RSD of 2
 replicates replicates Standard TE Curve: $y =$ Standard TE Curve: $y =$ Curve $348.31x + 1.979$ Curve
 $345.01x + 1.2758$ Equation TX Curve: $y =$ Equation TX Curve: $y = 279.85x + 6.8321$ $268.54x +$
 3.2308

TABLE-US-00046 TABLE A24 RiboGreen, Standard Curves, Post-Dialysis Samples Round 1
 Standard Set 1 Acceptance Criteria Measured Values Standard Curve R.sup.2 >0.975 >0.999
 Standard % RSD of 2 replicates $\leq 20\%$ <4.8% Standard Curve Equation TX Curve: $y = 293.84x -$
 2.242

TABLE-US-00047 TABLE A25 RiboGreen, Standard Curves, Post-Dialysis Samples Round 2
 Standard Set 1 Acceptance Criteria Measured Values Standard Curve R.sup.2 >0.975 >0.998
 Standard % RSD of 2 replicates $\leq 20\%$ <3.3% Standard Curve Equation TX Curve: $y = 267.84x +$
 1.4237

TABLE-US-00048 TABLE A26 RiboGreen, Standard Curves, Post-Dialysis Samples Round 3
 Plate 1 Acceptance Criteria Measured Values Standard Curve R.sup.2 >0.975 >0.996 Standard %
 RSD of 2 replicates $\leq 20\%$ <6.8% Standard Curve Equation TX Curve: $y = 214.45x + 11.38$

[0141] The stability of test sample was evaluated for three rounds of screening by visual
 appearance, dynamic light scattering (DLS), and RiboGreen RNA quantitation. Formulations
 containing 9% sucrose, 50 mM NaCl, and 50 mM Tris-HCl consistently displayed elegant,
 lyophilized cakes across all three rounds. Dextran 40 and iodixanol-containing formulations also
 produced elegant cakes with few minor bubbles. Conversely, formulations containing sodium
 benzoate or potassium sorbate without sucrose or iodixanol produced cakes which were collapsed
 or had many large bubbles. Upon reconstitution, sucrose-containing formulations displayed
 comparable to slightly increased opalescence in comparison to their pre-lyophilized states across
 three rounds. Particularly, following reconstitution with RNase-free water in Round 3, sucrose

formulations containing arginine and PVA10 were the least opalescent, whereas Dextran 40 formulations displayed the most opalescence.

[0142] Across three rounds, the DLS results following lyophilization revealed that formulations containing 9% sucrose and PVA10 showed smaller increases in hydrodynamic diameters compared to other formulations such as those containing iodixanol. In Rounds 1 and 2, formulations without P188 showed larger increases in hydrodynamic diameters after lyophilization than their P188-containing counterparts. In Round 3, Dextran 40 formulations showed larger increases in hydrodynamic diameters after lyophilization than other reconstituted formulations. Sucrose-containing formulations containing arginine and PVA10 were among those with the smallest increase in hydrodynamic diameters compared to their pre-lyophilization states, indicating these excipients optimally assist in the retention of test sample's hydrodynamic morphology.

[0143] Analysis of the RiboGreen RNA quantification at each round provided potential excipients for the next round of screening. In Round 2, reconstitution with MQ water resulted in slightly higher percent RNA encapsulation than reconstitution with TNS buffer. Furthermore, NaCl-containing formulations were found to provide improved encapsulation compared to formulations without NaCl in Round 2. Despite reported large increases in hydrodynamic diameters, iodixanol-containing formulations experienced smaller losses of encapsulation following lyophilization. In Round 3, Dextran 40 formulations displayed the largest loss of encapsulation and the greatest increase in hydrodynamic diameters. Formulations containing arginine and PVA10 showed high percentages of encapsulation after reconstitution. High concentrations of Tris-HCl also maintained high percentages of encapsulation yet showed large increases in hydrodynamic diameters after lyophilization. Therefore, arginine and PVA10 together with NaCl are better suited for higher RNA encapsulation following lyophilization compared to other formulations tested, maintaining their advantages to stability as previously observed by DLS and visual appearance. Formulation of 9% sucrose, 50 mM NaCl, 50 mM Tris-HCl, 100 mM arginine, and 0.05% PVA10 provides better stabilization for test sample through lyophilization and subsequent reconstitution by RNase-free water compared to other conditions tested.

[0144] In order to determine the improvements in excipient performance as described herein, stabilization during lyophilization was analyzed. Following dialysis, samples undergo spike with respective excipients targeting the various formulations, as shown in Table 21. The total volume of each formulation was targeted at 300 μ L. 100 μ L samples served as a pre-lyophilization control and analyzed without freezing. Two 96-well plates were prepared for lyophilization, with 100 μ L samples in each well.

[0145] A conservative lyophilization cycle was used to accommodate sufficient drying of all candidate formulations (Table 20). Following lyophilization, one lyophilized plate was sealed and stored at 5° C. for potential further testing. The other plate was visually inspected in the lyophilized state, and then reconstituted with 100 μ L of Milli-Q water. Both pre-lyophilization samples and reconstituted lyophilized samples were analyzed by visual inspection, DLS, and RiboGreen.

TABLE-US-00049 TABLE 20 Cycle Parameter Chamber Temperature Time Ramp Rate Pressure
Step (° C.) (min) (° C./min) (mT) Loading 5 N/A N/A N/A Freezing 5 to -50 55 1 N/A -50 120
N/A N/A -50 to -20 60 0.5 N/A -20 120 N/A N/A Primary Drying -20 1000* N/A 100
Secondary -20 to 20 80 0.5 100 Drying 20 600 N/A 100

TABLE-US-00050 TABLE 21 Excipient Screening Matrix Buffer Hydrophobic salts/AA Form #
(50 mM) Polymers (100 mM) NaCl Sucrose/Iodixanol 1 HEPES 0.05% Sodium Benzoate 50 mM
None 2 PVA10 Potassium Sorbate 3 0.2% Sodium Benzoate 4 P188 Potassium Sorbate 5 Tris-HCl
0.05% Potassium Sorbate 6 PVA10 Sodium Benzoate 7 0.2% Potassium Sorbate 8 P188 Sodium
Benzoate 9 HEPES 0.05% Potassium Sorbate None PVA10 10 0.2% Sodium Benzoate P188 11
Tris-HCl 0.05% Sodium Benzoate PVA10 12 0.2% Potassium Sorbate 13 P188 None 50 mM 2%
Sucrose 14 9% Sucrose 15 0.05% PVA10 16 0.2% Arginine None 17 P188 100 mM Arginine + 100
mM Potassium Sorbate 18 100 mM Proline + 100 mM Sodium Benzoate 19 HEPES 0.05% Proline

None 9% Iodixanol 20 PVA10 Arginine 21 0.2% Proline 9% Sucrose 22 P188 Arginine 23 Tris-HCl 0.05% Proline 9% Iodixanol 24 PVA10 Arginine 25 0.2% Proline 9% Sucrose 26 P188 Arginine 27 Proline 50 mM 2% Sucrose 28 20% Sucrose 29 Arginine 2% Sucrose 30 9% Sucrose [0146] Test sample was dialyzed against two base buffers (100 mM HEPES and 100 mM Tris, pH adjusted to 7.4) using Thermo Scientific™ Slide-A-Lyzer™, 20 kDa MWCO dialysis cassettes. At least a 10,000-fold buffer exchange was performed. Following dialysis, the pH of the dialyzed bulk DS is confirmed. RiboGreen analysis was performed to determine the Total RNA concentration of each formulation. Each dialyzed DP pool was diluted to 0.5 mg/mL with corresponding base buffers (100 mM HEPES or 100 mM Tris, at pH 7.4) and filtered. Excipient stock solutions were prepared in water and filtered, as shown in Table 22. The pH of the amino acid stock solutions is adjusted to pH 7.4. The 2× excipient solutions were prepared by combining all the necessary excipients for each formulation listed in Table 21. The 2× stock solutions are then spiked into the dialyzed formulations to reach the excipient parameters listed in Table 21 at 0.25 mg/mL test sample concentration.

TABLE-US-00051 TABLE 22 2 X Excipient Stock Solutions Volume in 1 mL 2x Reagent Stock Con. Target Con. Stock Sol'n Poloxamer 188 20% (w/v) 0.20% 20 µL PVA 10 4% (w/v) 0.05% 25 µL Sodium Benzoate 2M 0.1M 100 µL Potassium 2M 0.1M 100 µL Sorbate Sucrose 100% (w/v) 2%, 9%, 20% 40 µL, 180 µL, 400 µL Iodixanol 60% 9% 300 µL Proline 1M 0.1M 200 µL Arginine 1M 0.1M 200 µL NaCl 2.5M 0.05M 40 µL

[0147] A preparation was prepared using 100 µL samples of each formulation which served as a pre-lyophilization control and was analyzed without freezing. Additionally, two 96-well plates were prepared for lyophilization, with 100 µL samples in each well. A conservative lyophilization cycle was used to accommodate sufficient drying of all candidate formulations. Following lyophilization, one lyophilized plate was sealed and stored at 5° C. for further testing. The other plate was visually inspected in the lyophilized state, and then reconstituted with 100 µL of Milli-Q water. Both pre-lyophilization samples and reconstituted lyophilized samples were analyzed by visual inspection, DLS, and RiboGreen.

[0148] Excipient screening during Round 2 was performed to assess excipients that can stabilize test sample during lyophilization. Following dialysis, samples may undergo a spike with respective excipients targeting the various formulations, as shown in Table 24. Thirty (30) formulations were prepared, targeting a 0.25 mg/mL final concentration of test sample. The total volume of each formulation was 400 µL. 100 µL samples serve as a pre-lyophilization control and were analyzed without freezing. Of note, before dialysis, the drug product was aliquoted and kept at 5° C. without refreezing. The DP control was analyzed together with the pre-lyophilization controls. Three 96-well plates were prepared for lyophilization, with 100 µL samples in each well. A revised lyophilization cycle is used to accommodate sufficient drying of all candidate formulations (Table 23). Following lyophilization, one lyophilized plate was sealed and stored at 5° C. for potential further testing. The other two plates were visually inspected in the lyophilized state, and then reconstituted with 100 µL of Milli-Q water and 100 µL of 12.5 mM Tris, 12.5 mM NaCl, 2.25% sucrose, pH 7.5, respectively. Both pre-lyophilization samples and reconstituted lyophilization samples were analyzed by visual inspection, DLS, and RiboGreen.

TABLE-US-00052 TABLE 23 Cycle Parameter Chamber Temperature Time Ramp Rate Pressure Step (° C.) (min) (° C./min) (mT) Loading 5 N/A N/A N/A Freezing 5 to -50 55 1 N/A -50 120 N/A N/A -50 to -30 40 0.5 N/A -30 120 N/A N/A Primary Drying -30 1000* N/A 100 Secondary -30 to 5 70 0.5 100 Drying 5 1200 N/A 100

TABLE-US-00053 TABLE 24 Excipient Screening Matrix Sucrose/ Buffer PVA10 P188 Amino Amino Acid Iodixanol NaCl Form # (50 mM) (%) (%) Acid Con. (M) (9%) (mM) 1 Tris-HCl 0.1 0.4 Arginine 0.05 Sucrose NA 2 0.1 0.2 Proline 0.05 Sucrose NA 3 0.05 0.4 Proline 0.1 Sucrose NA 4 0.05 NA Arginine 0.1 Sucrose NA 5 0.035 0.2 Proline 0.1 Sucrose NA 6 0.035 NA Proline 0.1 Sucrose NA 7 0.1 0.2 Arginine 0.05 Iodixanol NA 8 0.1 NA Proline 0.1 Iodixanol NA 9 0.05 0.4

Arginine 0.1 Iodixanol NA 10 0.05 0.2 Proline 0.05 Iodixanol NA 11 0.035 0.4 Proline 0.05 Iodixanol NA 12 0.035 NA Arginine 0.05 Iodixanol NA 13 0.1 0.4 Arginine 0.05 Sucrose 50 14 0.1 0.2 Proline 0.05 Sucrose 50 15 0.05 0.4 Proline 0.1 Sucrose 50 16 0.05 NA Arginine 0.1 Sucrose 50 17 0.035 0.2 Proline 0.1 Sucrose 50 18 0.035 NA Proline 0.1 Sucrose 50 19 0.1 0.2 Arginine 0.05 Iodixanol 50 20 0.1 NA Proline 0.1 Iodixanol 50 21 0.05 0.4 Arginine 0.1 Iodixanol 50 22 0.05 0.2 Proline 0.05 Iodixanol 50 23 0.035 0.4 Proline 0.05 Iodixanol 50 24 0.035 NA Arginine 0.05 Iodixanol 50 25 NA 0.2 Proline 0.1 Iodixanol 50 26 NA 0.2 Arginine 0.1 Iodixanol 50 27 NA NA Proline 0.1 Iodixanol 50 28 NA NA Arginine 0.1 Sucrose 50 29 NA NA NA NA Iodixanol 50 30 NA NA NA NA NA Sucrose 50

[0149] The Experimental procedure included in Round 3 excipient studies were performed using visual, DLS and RiboGreen. Formulation preparation of test sample was dialyzed against base buffers (100 mM Tris, pH adjusted to 7.4) using Thermo Scientific™ Slide-A-Lyzer™, 20 kDa MWCO dialysis cassettes. At least a 10,000-fold buffer exchange was performed. Following dialysis, the pH of the dialyzed bulk DS was confirmed. RiboGreen analysis was performed to determine the Total RNA concentration of each formulation. The dialyzed DP pool was diluted to 0.5 mg/mL with the base buffer and filtered. Excipient stock solutions were prepared in water and filtered, as shown in Table 25. The pH of the amino acid stock solutions was adjusted to a pH 7.4. The 2× excipient solutions were prepared by combining all the necessary excipients for each formulation listed in Table 24. The 2× stock solutions are then spiked into the dialyzed formulations to reach the excipient parameters listed in Table 24 at 0.25 mg/mL test sample concentration.

TABLE-US-00054 TABLE 25 2 X Excipient Stock Solutions Volume in 1 mL 2x Reagent Stock Con. Target Con. Stock Solution Poloxamer 188 20% (w/v) 0.2%, 0.4% 20 µL, 40 µL PVA 10 4% (w/v) 0.035%, 0.05%, 0.1% 17.5 µL, 25 µL, 50 µL Sucrose 60% (w/v) 9% 300 µL Iodixanol 60% 9% 300 µL Proline 1M 0.05M, 0.1M 100 µL, 200 µL Arginine 1M 0.05M, 0.1M 100 µL, 200 µL NaCl 2.5M 0.05M 40 µL

[0150] Following preparation, 100 µL samples of each formulation served as a pre-lyophilization control and were analyzed without freezing. Three 96-well plates were prepared for lyophilization, with 100 µL samples in each well. Following lyophilization, one lyophilized plate was sealed and stored at 5° C. for potential further testing. The other two plates were visually inspected in the lyophilized state, and then reconstituted with 100 µL of two reconstitution buffers, respectively. The two reconstitution buffers are 1) Milli-Q water and 2) 12.5 mM Tris, 12.5 mM NaCl, 2.25% sucrose, pH 7.5. Both pre-lyophilization samples and reconstituted lyophilization samples were analyzed by visual inspection, DLS, and RiboGreen.

[0151] Round 3 included the respective excipients targeting the various formulations desired to be utilized for therapeutic formulations within the scope of the described invention. As shown in Table 27, thirty (30) formulations were prepared, targeting a 0.25 mg/mL final concentration of test sample. The total volume of each formulation was 400 µL. 100 µL samples served as a pre-lyophilization control and were analyzed without freezing. Of note, before dialysis, the drug product was aliquoted and kept at 5° C. without refreezing. The DP control was analyzed together with the pre-lyophilization controls. Three 96-well plates were prepared for lyophilization, with 100 µL samples in each well. A revised lyophilization cycle was used to accommodate sufficient drying of all candidate formulations (Table 26).

[0152] Following lyophilization, two lyophilized plates were sealed and stored at 5° C. The other plate was visually inspected in the lyophilized state, and then reconstituted with 100 µL of Milli-Q water. Both pre-lyophilization samples and reconstituted lyophilization samples were analyzed by visual inspection, DLS, and RiboGreen.

TABLE-US-00055 TABLE 26 Cycle Parameter Temperature Time Ramp Rate Chamber Step (° C.) (min) (° C./min) Pressure (mT) Loading 5 N/A N/A N/A Freezing 5 to -50 55 1 N/A -50 120 N/A N/A -50 to -30 40 0.5 N/A -30 120 N/A N/A Primary Drying -30 1000* N/A 100 Secondary

chloride, 0.1% polyvinyl alcohol, 15% sucrose at pH 7.6. Bulk product was sterile filtered through a 0.22 µm PES membrane.

[0158] 1 mL aliquots of Lots 230607A and B were filled into 3 mL glass vials and partially stoppered with West Pharmaceutical lyo stoppers. Vials were lyophilized using the cycle described in Table 29.

[0159] 0.5 mL aliquots of Lots 230523A and B were filled into 3 mL glass vials and partially stoppered with West Pharmaceutical lyo stoppers. Vials were lyophilized using the cycle described in Table 30.

TABLE-US-00058 TABLE 29 Lyophilization Cycle Parameters Used for Lots 230607A and B

Chamber Step	Temperature (° C.)	Time	Ramp Rate	Pressure	Loading
NA	18° C.	10 min	5° C./min	NA	NA
Freezing	18° C. to 5° C.	13 min	1° C./min	NA	5° C.
NA	5° C.	60 min	NA	NA	5 to -50° C.
Maximum NA	NA	-50° C.	300 min	NA	NA
Annealing	-50 to -10° C.	40 min	1° C./min	NA	-10° C.
NA	-10° C.	180 min	NA	NA	-10° C. to -50° C.
Primary	-50° C.	40 min	1° C./min	NA	-50° C.
NA	-50° C.	180 min	NA	NA	-50° C. to -30° C.
Drying	-30° C.	200 min	0.1° C./min	100 mTorr	NA
Secondary	-30° C.	2400 min	NA	100 mTorr	-30° C. to -24° C.
Drying	-24° C.	60 min	0.1° C./min	100 mTorr	-24° C.
Secondary	-24° C.	900 min	NA	100 mTorr	-24° C. to 25° C.
Drying	25° C.	490 min	0.1° C./min	100 mTorr	NA

TABLE-US-00059 TABLE 30 Lyophilization Cycle Parameters Used for Lots 230523A and B

Chamber Step	Temperature (° C.)	Time	Ramp Rate	Pressure	Loading
NA	18° C.	NA	NA	NA	NA
Freezing	18° C. to 5° C.	13 min	1° C./min	NA	5° C.
NA	5° C.	30 min	NA	NA	5 to -50° C.
Maximum NA	NA	-50° C.	55 min	1° C./min	NA
Annealing	-50 to -15° C.	35 min	1° C./min	NA	-15° C.
NA	-15° C.	180 min	NA	NA	-15° C. to -50° C.
Primary	-50° C.	35 min	1° C./min	NA	-50° C.
NA	-50° C.	180 min	NA	NA	-50° C. to -30° C.
Drying	-30° C.	200 min	0.1° C./min	100 mTorr	NA
Secondary	-30° C.	2400 min	NA	100 mTorr	-30° C. to -22° C.
Drying	-22° C.	80 min	0.1° C./min	100 mTorr	-22° C.
Secondary	-22° C.	600 min	NA	100 mTorr	-22° C. to 25° C.
Drying	25° C.	470 min	0.1° C./min	100 mTorr	NA

[0160] Following lyophilization, samples were placed on stability at -20° C., 5° C., 25° C., and 40° C.

[0161] Lots 230607A and B were reconstituted and analyzed by DLS, encapsulation efficiency, and mRNA purity by CE-LIF at the time points listed in Table 31. Potency was analyzed at TO and 2 week at 40° C., 4 week 25° C., 6 months at -20° C., and 5° C. for lot 230607A and 1 and 6 months at -20° C. and 5° C. for lot 230607B (Table 35).

[0162] For the potency method, HepG2 AGL knockout (KO) cells are seeded into a 96-well plate and cultured overnight. On the same day, the mRNA-LNP test samples, reference standard, and assay control are diluted to the target concentration and then incubated with human serum to allow desorption of the PEG lipid from the LNP surface. The next day, serial dilutions of the human serum-treated DP samples are prepared and then incubated with the HepG2 AGL KO cells overnight to allow LNP entry into the cells and translation of the GDE protein from the delivered mRNA. On the third day, the cells are fixed, permeabilized, blocked, and incubated with a specific anti-GDE antibody. On the last day, a secondary antibody labeled with a fluorescent tag (800 nm) along with a CellTag stain (700 nm) used to normalize the cell numbers are added to stain the cells. The final stained 96-well plate is scanned on the LI-COR Odyssey CLx imaging system at 800 nm and 700 nm. The normalized ICW signals are calculated by dividing the GDE expression (800 nm) by the total cell number (700 nm). Finally, dose-response curves are generated for the test samples, reference standard, and control sample. Relative potency is calculated by comparing the EC50 of the test samples to the reference standard using a sigmoidal 4PL curve fitting model in the PLA software.

[0163] Lots 230523A and B were analyzed by DLS, encapsulation efficiency, mRNA purity by CE-LIF at the time point listed in Table 32. Potency was analyzed at TO and 1 week at 40° C., 6 months at -20° C., 5° C. and 25° C. (Table 35).

TABLE-US-00060 TABLE 31 Stability testing plan for lots 230607A and B

Purity by CE-LIF	DLS	40° C.	25° C.	5° C.	-20° C.	40° C.	25° C.	5° C.	-20° C.	Encapsulation Efficiency	T0	x	x	40° C.
------------------	-----	--------	--------	-------	---------	--------	--------	-------	---------	--------------------------	----	---	---	--------

25° C. 5° C. -20° C. 1 day x x x x x 2 days x x x x x 1 week x x x x x 2 weeks x x x x 4 weeks x x 3 months x x 6 months x x x x x .sup.1 Lot 230607B 5° C. 2 week time point was not tested for Purity by CE-LIF

TABLE-US-00061 TABLE 32 Stability testing plan for lots 230523A and B 40° C. 25° C. 5° C. -20° C. T0 x 1 week x x 2 weeks x 3 months x 4 months x x 6 months x x x

[0164] Lots 230523A and B remained stable by polydispersity index [PDI] (FIG. 20) and encapsulation efficiency (FIG. 22) for 6 months at 5° C., 6 months -20° C., 6 months at 25° C. and 2 weeks at 40° C. There is a slight increase in size observed for lots stored 1 week at 40° C. (FIG. 19) but all other conditions remained stable. FIG. 24 shows that there's an increasing loss of purity with increasing temperature storage with all conditions showing a loss of purity, but it appears to be independent of NP ratio.

[0165] Lots 230607A and B size remained stable after 1 week at 40° C., 1 weeks at 25° C., and 6 months at 5° C. and -20° C. (FIG. 21 and Table 33) and by PDI after 2 weeks at 40° C., 1 month at 25° C., 6 months at 5° C., and 6 months at -20° C. (Table 34). Purity decreased under accelerate and stress storage conditions, 25° C. and 40° C. respectively, with a higher rate of degradation observed at 40° C. Purity remained unchanged after 6 months at 5° C. and -20° C. (FIG. 23).

[0166] It appears that the stability of the mRNA-LNP particle is dependent on the lyophilization cycle as lots 230607A and B have a slower rate of degradation compared to lots 230523A and B. TABLE-US-00062 TABLE 33 Change in size from T0 for lots 230607 A (NP9) B (NP6) [nm] [nm]

6 M -20 C.	6 -2	6 M 5 C.	12	3	1 M 25 C.	10	11	2 W 40 C.	16	17
------------	------	----------	----	---	-----------	----	----	-----------	----	----

TABLE-US-00063 TABLE 34 Change in PDI from T0 for lots 230607 A (NP9) B (NP6) 6 M -20 C. -0.01 0.01 6 M 5 C. 0.02 0.01 1 M 25 C. 0.00 0.00 2 W 40 C. -0.03 -0.01

TABLE-US-00064 TABLE 35 Potency from T0 for lots 230607A and B and 230523A and B -20 C. -20 C. 5 C. 5 C. 25 C. 25 C. 40 C. 40 C. 0 1 M 6 M 1 M 6 M 28 D 6 M 1 W 2 W

230607A(NP9)	137	267	408	206	391	115	13	230607B (NP6)	133	287	306	230523A (NP9)	96	92
51	0	0	230523B (NP6)	72	86	49	0	0						

ENUMERATED EMBODIMENTS

[0167] Although the present disclosure has been described with respect to one or more particular embodiments, it will be understood that other embodiments of the present disclosure may be made without departing from the scope of the present disclosure.

[0168] Further examples and embodiments of the present disclosure are provided in the enumerated clauses which follow:

[0169] 1. A composition comprising a cargo molecule and lipid nanoparticle, wherein the composition is made using the method of any one of clauses 22-88.

[0170] 2. The composition of clause 1, wherein the cargo molecule comprises a nucleic acid.

[0171] 3. The composition of clause 2, wherein the nucleic acid comprises a DNA or an RNA, or an antisense oligonucleotide (ASO).

[0172] 4. The composition of clause 3, wherein the nucleic acid comprises a miRNA, mRNA, siRNA, ASO, SAM, circular RNA, or a combination thereof.

[0173] 5. The composition of clause 3, wherein the RNA is mRNA.

[0174] 6. The composition of any one of clauses 2-5, wherein the nucleic acid is from about 14 nucleotides to about 15,000 nucleotides in length.

[0175] 7. The composition of any one of clauses 1-6, wherein a total lipid to nucleic acid weight ratio in the composition is about 50:1 to about 4:1, or 40:1 to about 25:1, or 35:1 to about 20:1, or 20:1 to about 10:1, or 15:1 to about 4:1.

[0176] 8. The composition of any one of clauses 1-7, wherein the lipid nanoparticle in the composition after being lyophilized and then reconstituted is increased in size between 3-18% relative to the size of the lipid nanoparticle prior to lyophilization and reconstitution.

[0177] 9. The composition of any one of clauses 1-7, wherein the lipid nanoparticle in the composition after being lyophilized and then reconstituted is increased in size no more than 17%

relative to the size of the lipid nanoparticle prior to lyophilization and reconstitution.

[0178] 10. The composition of any one of clauses 2-9, wherein the composition after lyophilization comprises a nucleic acid encapsulation efficiency (EE) of at least 70%.

[0179] 11. The composition of any one of clauses 2-9, wherein the composition after lyophilization comprises a nucleic acid encapsulation efficiency (EE) of at least 80%.

[0180] 12. The composition of any one of clauses 2-9, wherein the composition after lyophilization comprises a nucleic acid encapsulation efficiency (EE) of at least 85%.

[0181] 13. The composition of any one of clauses 2-9, wherein the composition after lyophilization comprises a nucleic acid encapsulation efficiency (EE) of at least 90%.

[0182] 14. The composition of any one of clauses 1-13, wherein the cargo molecule comprises mRNA with a concentration of at least 0.05 mg/ml.

[0183] 15. The composition of clause 14, wherein the mRNA is combined with a dilution buffer comprising Ala, Arg, His, Asp, Pro, Glu, Met, Thr, or any combination thereof.

[0184] 16. The composition of clause 15, wherein the dilution buffer comprises at least 2 amino acids selected from a group comprising Ala, Arg, His, Asp, Pro, Glu, Met, Thr, or a combo thereof, wherein one of the amino acids is positively charged and one amino acid is negatively charged.

[0185] 17. The composition of any one of clauses 1-16, wherein the composition comprises an excipient in a weight ratio of excipient to cargo molecule of about 25:1 to 400:1.

[0186] 18. The composition of any one of clauses 1-17, wherein the excipient is in a weight ratio of excipient to cargo molecule of about 2:1 to about 10:1.

[0187] 19. The composition of any one of clauses 1-17, wherein the excipient is in a weight ratio of excipient to cargo molecule of about 3:1 to about 8:1.

[0188] 20. The composition of any one of clauses 1-17, wherein the excipient is in a weight ratio of excipient to cargo molecule of about 1:1 to about 12:1.

[0189] 21. The composition of any one of clauses 1-20, wherein the excipient is selected from: citrate, tris, arginine, NaCl, sucrose, PVA, phosphate, HEPES, trehalose, KCl, acetate, bis-tris, histidine, glucose, lactose, raffinose, alanine, asparagine, proline, glutamic acid, methionine, threonine, and combinations thereof.

[0190] 22. A method of making a composition comprising a lipid nanoparticle and a cargo molecule disposed therein, wherein the composition is capable of being lyophilized, the method comprising: [0191] a. combining an aqueous solution of the cargo molecule and a lipid composition with an ionizable organic solvent solution, thereby forming a lipid/cargo molecule suspension.

[0192] 23. The method of clause 22, wherein the aqueous solution and/or lipid composition comprises an organic solvent.

[0193] 24. The method of any one of clauses 22-23, wherein the cargo molecule comprises a nucleic acid.

[0194] 25. The method of any one of clauses 22-24, wherein the nucleic acid comprises a DNA or an RNA, or an ASO.

[0195] 26. The method of clause 25, wherein the nucleic acid comprises a miRNA, mRNA, SiRNA, ASO, SAM, circular RNA, or a combination thereof.

[0196] 27. The method of any one of clauses 25-26, wherein the RNA is mRNA.

[0197] 28. The method of any one of clauses 22-27, wherein the lipid composition comprises an ionizable lipid composition.

[0198] 29. The method of any one of clauses 22-28, wherein the organic solvent comprises at least one of DSPC, cholesterol, PEG-DMG, an ionizable lipid, and combinations thereof.

[0199] 30. The method of any one of clauses 22-29, further comprising the step of adding an excipient to the aqueous solution of the cargo molecule prior to being combined with the lipid composition.

[0200] 31. The method of any one of clauses 22-29, further comprising the step of adding an

excipient to the lipid/cargo molecule suspension.

[0201] 32. The method of any one of clauses 30-31, wherein the excipient is selected from citrate, tris, arginine, NaCl, sucrose, PVA, phosphate, HEPES, trehalose, KCl, acetate, bis-tris, histidine, glucose, lactose, raffinose, alanine, asparagine, proline, glutamic acid, methionine, threonine, and combinations thereof.

[0202] 33. The method of any one of clauses 22-32, wherein the ionizable organic solvent solution comprises at least one of an ionizable lipid, DSPC, cholesterol, PEG-DMG and combinations thereof.

[0203] 34. The method of clause 33, wherein the ionizable organic solvent solution comprises a N/P molar ratio ≥ 2.5 .

[0204] 35. The method of clause 33, wherein ionizable organic solvent solution comprises a N/P molar ratio ≥ 6.0 .

[0205] 36. The method of clause 33, wherein the ionizable organic solvent solution comprises a N/P molar ratio ≥ 9.0 .

[0206] 37. The method of any one of clauses 23-36, further comprising the step of (b) diluting the lipid/cargo molecule suspension combined with the organic solvent to lower the concentration of organic solvent to $\leq 10\%$ of the total composition.

[0207] 38. The method of any one of clauses 23-36, further comprising the step of (b) diluting the lipid/cargo molecule suspension combined with the organic solvent to lower the concentration of organic solvent to $\leq 3\%$ of the total composition.

[0208] 39. The method of any one of clauses 22-38, further comprising the step of increasing the pH of the composition to a pH at or above 5.0.

[0209] 40. The method of clause 39, further comprising the step (b) of clause 37 or 38, wherein the pH is increased after step (b).

[0210] 41. The method of clause 39, further comprising the step (b) of clause 37 or 38, wherein the pH is increased prior to step (b).

[0211] 42. The method of any one of clauses 40 or 41, wherein the pH is at least 6.

[0212] 43. The method of any one of clauses 40 or 41, wherein the pH is between 6.9 and 7.9.

[0213] 44. A method of lyophilizing a composition comprising the lipid/cargo molecule suspension of any one of clauses 22-43.

[0214] 45. The method of clause 44, wherein the lyophilized composition comprises the lipid nanoparticle having the cargo molecule disposed therein, thereby forming a lyophilized cargo molecule/LNP particle, the method comprising freezing the lipid/cargo molecule suspension.

[0215] 46. The method of any one of clauses 44 or 45, further comprising the step of drying the lipid/cargo molecule suspension, wherein the drying step is performed by a process selected from a group comprising spin-freeze drying, continuous freeze-drying, sublimation, desorption, vacuum foam drying, spray drying, evaporation and combination thereof.

[0216] 47. The method of clause 46, wherein the drying step comprises sublimation to achieve a moisture of $\leq 5\%$.

[0217] 48. The method of any one of clauses 44-47, wherein the potency and/or purity yield of the lyophilized composition is $\geq 70\%$.

[0218] 49. The method of any one of clauses 44-47, wherein the potency and/or purity yield of the lyophilized composition is $\geq 80\%$.

[0219] 50. The method of any one of clauses 44-47, wherein the potency and/or purity yield of the lyophilized composition is $\geq 90\%$.

[0220] 51. The method of any one of clauses 44-47, wherein the potency and/or purity yield of the lyophilized composition is $\geq 100\%$.

[0221] 52. The method of any one of clauses 44-51, wherein the lyophilized cargo molecule/LNP particle size is < 100 nm.

[0222] 53. The method of any one of clauses 44-51, wherein the lyophilized cargo molecule/LNP

particle size is <90 nm.

[0223] 54. The method of any one of clauses 44-51, wherein the lyophilized cargo molecule/LNP particle size is <80 nm.

[0224] 55. The method of any one of clauses 44-51, wherein the lyophilized cargo molecule/LNP particle size is <75 nm.

[0225] 56. The method of any one of clauses 44-55, wherein the lyophilized cargo molecule/LNP encapsulation efficiency (EE) is $\geq 80\%$.

[0226] 57. The method of any one of clauses 44-55, wherein the lyophilized cargo molecule/LNP encapsulation efficiency (EE) is $\geq 90\%$.

[0227] 58. The method of any one of clause 44-57, wherein the cargo molecule comprises a therapeutic agent at a starting concentration of ≤ 2 mg/ml under acidic conditions with a pH under 5, preferably between 3-4.5.

[0228] 59. The method of any one of clause 44-57, wherein the cargo molecule is an mRNA at a starting concentration of ≤ 1 mg/ml under acidic conditions with a pH under 5, preferably between 3-4.5.

[0229] 60. The method of any one of clause 44-57, wherein the cargo molecule is an mRNA at a starting concentration of ≤ 0.55 mg/ml under acidic conditions with a pH under 5, preferably between 3 and 4.5.

[0230] 61 The method of any one of clause 44-57, wherein the cargo molecule is an mRNA at a starting concentration of ≤ 0.22 mg/ml under acidic conditions with a pH under 5, preferably between 3 to 4.5.

[0231] 62. The method of any one of clauses 22-61, wherein the cargo molecule, in a non-lyophilized state, prior to encapsulation, is dissolved in a buffer under acidic pH conditions below the pKa of the lipid composition and wherein the buffer optionally comprises 50 mM NaCl.

[0232] 63. The method of any one of clauses 22-61, wherein the cargo molecule, in a non-lyophilized state, prior to encapsulation, is dissolved in a buffer under acidic pH conditions below the pKa of the lipid composition and wherein the buffer comprises a salt between 25 mM and 100 mM, wherein the salt is preferably NaCl.

[0233] 64 The method of any one of clauses 30-63, wherein the excipient(s) comprises at least one positively charged and at least one negatively charged compound.

[0234] 65. The method of clause 64, wherein the excipient comprise at least two amino acids selected from Ala, Arg, His, Asp, Pro, Glu, Met, Thr, and combinations thereof.

[0235] 66. The method of any one of clauses 38-65, wherein the dilution step comprises a first dilution step to increase the pH to about 6 and a second dilution step to increase the pH to between 6.8 and 7.8.

[0236] 67. The method of clause 66, wherein the first dilution step comprises at least one of phosphate, HEPES, tris, bis-tris, acetate, citrate, NaCl, KCl, glutamic acid, arginine, histidine, methionine, glucose, sucrose, lactose, trehalose, raffinose, and combinations thereof.

[0237] 68. The method of any one of clauses 66 or 67, wherein the second dilution buffer comprises at least one buffer compound selected from amino acids, sugars, salts and combinations thereof.

[0238] 69. The method of any one of clauses 22-68, further comprising the step of adding a stabilizer to the aqueous solution of the cargo molecule.

[0239] 70. The method of clause 69, wherein the stabilizer is added prior to mixing the cargo molecule with the organic solvent and wherein the stabilizer is optionally present in any one of the at least one dilution steps.

[0240] 71. The method of any one of clauses 30-70, wherein the excipient comprises a hydrophilic polymeric compound.

[0241] 72. The method of any one of clauses 30-71, wherein the excipient further comprises a plasticizer.

[0242] 73. A method of preparing a reconstituted composition comprising a lipid nanoparticle having a cargo molecule disposed therein from a lyophilized composition prepared according to the method of any one of clauses 44-72, comprising the step of reconstituting the lyophilized composition with a diluent.

[0243] 74. The method of clause 73, wherein the lipid nanoparticle in the composition after being lyophilized and then reconstituted is increased in size between 3-18% or to a size $\leq 18\%$ relative to the size of the lipid nanoparticle prior to lyophilization.

[0244] 75. The method of clause 73, wherein the lipid nanoparticle in the composition after being lyophilized and then reconstituted is increased in size no more than 17% relative to the size of the lipid nanoparticle prior to lyophilization.

[0245] 76. The method of any one of clauses 73-75, wherein the cargo molecule comprises an RNA and the encapsulation efficiency (EE) of the RNA is at least 70%.

[0246] 77. The method of any one of clauses 73-75, wherein the cargo molecule comprises an RNA and the encapsulation efficiency (EE) of the RNA is at least 80%.

[0247] 78. The method of any one of clauses 73-75, wherein the cargo molecule comprises an RNA and the encapsulation efficiency (EE) of the RNA is at least 85%.

[0248] 79. The method of any one of clauses 73-75, wherein the cargo molecule comprises an RNA and the encapsulation efficiency (EE) of the RNA is at least 90%.

[0249] 80. The method of any one of clauses 73-79, wherein the encapsulation efficiency (EE) of the composition after being reconstituted has an average variance of no more than 20% of the lyophilized or aqueous form of the composition.

[0250] 81. The method of any one of clauses 73-79, wherein the encapsulation efficiency (EE) of the composition after being reconstituted has an average variance of between 10-20% of the lyophilized or aqueous form of the composition.

[0251] 82. A method of preparing a lyophilized composition comprising lipid nanoparticles with reduced percentage of water in a lipid nanoparticle core, the method comprising: [0252] (a) combining an aqueous solution of a cargo molecule with a lipid composition, thereby forming a lipid/cargo molecule suspension; [0253] (b) mixing the lipid/cargo molecule suspension with an organic solvent; [0254] (c) adding an excipient, wherein the excipient is added before or after the aqueous solution of the cargo molecule is combined with the lipid composition; [0255] (d) diluting the lipid/cargo molecule suspension combined with the organic solvent to lower the concentration of organic solvent to $\leq 10\%$ of the total composition by volume; thereby forming a lipid nanoparticle composition comprising lipid nanoparticles comprising a shell of the lipid composition and a core comprising the cargo molecule in aqueous solution with the organic solvent; [0256] (e) lyophilizing the lipid nanoparticle composition; and [0257] (f) drying the lipid nanoparticle composition, thereby forming the lyophilized composition with lipid nanoparticles having reduced percentage of water in the lipid nanoparticle core.

[0258] 83. The method of clause 82, wherein the moisture content of the lyophilized composition is $\leq 5\%$ by weight.

[0259] 84. The method of clause 82 or 83, wherein the drying step is performed by a process selected from spin-freeze drying, continuous freeze-drying, sublimation, desorption, vacuum foam drying, spray drying, evaporation, and combinations thereof.

[0260] 85. A method of preparing a lipid nanoparticle composition comprising a lipid nanoparticle prepared according to the method of any one of clauses 82-84, the method comprising reconstituting the lyophilized composition in a diluent.

[0261] 86. The method clause 85, wherein the diluent comprises saline, water, 50-500 mM NaCl, arginine, D5W, sucrose, polyvinyl alcohol (PVA), tris, or combinations thereof.

[0262] 87. The method of any one of clauses 82-86, wherein the drying step occurs at a temperature between 30°C . to -50°C . 88. The method of any one of clauses 82-87, wherein the drying step is performed at a time range of between 10-3000 minutes.

[0263] 89. The composition of any one of clauses 1-21, wherein the cargo molecule comprises an mRNA encoding amylo-alpha-1, 6-glucosidase, 4-alpha-glucanotransferase (AGL), or a fragment thereof.

[0264] 90. The composition of clause 89, wherein the composition is for ameliorating, preventing, delaying onset, or treating a disease or condition associated with reduced activity of AGL in a subject.

[0265] 91. The composition of any one of clauses 1-21, wherein the cargo molecule comprises an mRNA encoding bile salt export pump (BSEP).

[0266] 92. The composition of clause 91, wherein the composition is for ameliorating, preventing, delaying onset, or treating a disease or condition associated with reduced activity of BSEP in a subject.

[0267] The above specification, examples and data provide a complete description of the manufacture and use of the methods and compositions of the invention. Since many embodiments of the invention can be made without departing from the spirit and scope of the invention, the invention resides in the claims hereinafter appended.

Claims

1. A method of making a composition comprising a lipid nanoparticle and a cargo molecule disposed therein, wherein the composition is capable of being lyophilized, the method comprising: (a) combining an aqueous solution of the cargo molecule and a lipid composition with an ionizable organic solvent solution, thereby forming a lipid/cargo molecule suspension.
2. The method of claim 1, wherein the aqueous solution and/or lipid composition comprises an organic solvent.
3. The method of claim 1, wherein the cargo molecule comprises a nucleic acid.
4. The method of claim 1, wherein the lipid composition comprises an ionizable lipid composition.
5. The method of claim 1, further comprising the step of adding an excipient to the aqueous solution of the cargo molecule prior to being combined with the lipid composition.
6. The method of claim 1, further comprising the step of adding an excipient to the lipid/cargo molecule suspension.
7. The method of claim 1, wherein the ionizable organic solvent solution comprises a N/P molar ratio ≥ 2.5 .
8. The method of claim 2, further comprising the step of (b) diluting the lipid/cargo molecule suspension combined with the organic solvent to lower the concentration of organic solvent to $\leq 10\%$ of the total composition.
9. The method of claim 1, further comprising the step of increasing the pH of the composition to a pH at or above 5.0.
10. The method of claim 1, wherein the cargo molecule, in a non-lyophilized state, prior to encapsulation, is dissolved in a buffer under acidic pH conditions below the pKa of the lipid composition.
11. A method of lyophilizing a composition comprising the lipid/cargo molecule suspension of claim 1, wherein the lyophilized composition comprises the lipid nanoparticle having the cargo molecule disposed therein, thereby forming a lyophilized cargo molecule/LNP particle, the method comprising freezing the lipid/cargo molecule suspension.
12. The method of claim 11, further comprising the step of drying the lipid/cargo molecule suspension, wherein the drying step is performed by a process selected from a group comprising spin-freeze drying, continuous freeze-drying, sublimation, desorption, vacuum foam drying, spray drying, evaporation and combination thereof.
13. The method of claim 12, wherein the lyophilized cargo molecule/LNP particle size is < 100 nm.
14. The method of claim 12, wherein the cargo molecule comprises a therapeutic agent at a starting

concentration of ≤ 2 mg/ml under acidic conditions with a pH under 5, preferably between 3-4.5.

15. A method of preparing a lyophilized composition comprising lipid nanoparticles with reduced percentage of water in a lipid nanoparticle core, the method comprising: (a) combining an aqueous solution of a cargo molecule with a lipid composition, thereby forming a lipid/cargo molecule suspension; (b) mixing the lipid/cargo molecule suspension with an organic solvent; (c) adding an excipient, wherein the excipient is added before or after the aqueous solution of the cargo molecule is combined with the lipid composition; (d) diluting the lipid/cargo molecule suspension combined with the organic solvent to lower the concentration of organic solvent to $\leq 10\%$ of the total composition by volume; thereby forming a lipid nanoparticle composition comprising lipid nanoparticles comprising a shell of the lipid composition and a core comprising the cargo molecule in aqueous solution with the organic solvent; (e) lyophilizing the lipid nanoparticle composition; and (f) drying the lipid nanoparticle composition, thereby forming the lyophilized composition with lipid nanoparticles having reduced percentage of water in the lipid nanoparticle core.

16. The method of claim 15, wherein: (a) the moisture content of the lyophilized composition is $\leq 5\%$ by weight; (b) the drying step occurs at a temperature between 30°C. to -50°C. ; and/or (c) the drying step is performed at a time range of between 10-3000 minutes.

17. A composition comprising a cargo molecule and lipid nanoparticle, wherein the composition is made using the method of claim 1.

18. The composition of claim 17, wherein the cargo molecule comprises a nucleic acid.

19. The composition of claim 18, wherein a total lipid to nucleic acid weight ratio in the composition is about 50:1 to about 4:1, or 40:1 to about 25:1, or 35:1 to about 20:1, or 20:1 to about 10:1, or 15:1 to about 4:1.

20. The composition of claim 19, wherein the cargo molecule comprises mRNA with a concentration of at least 0.05 mg/ml.

21. The composition of claim 20, wherein the mRNA is combined with a dilution buffer comprising Ala, Arg, His, Asp, Pro, Glu, Met, Thr, or any combination thereof.

22. The composition of claim 17, wherein the composition comprises an excipient in a weight ratio of excipient to cargo molecule of about 25:1 to 400:1.

23. The composition of claim 17, wherein the cargo molecule comprises an mRNA encoding amylo-alpha-1, 6-glucosidase, 4-alpha-glucanotransferase (AGL), or a fragment thereof.

24. The composition of claim 17, wherein the cargo molecule comprises an mRNA encoding bile salt export pump (BSEP).
